

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
3. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.
4. Social Workers: For additional information regarding standards and indicators for cultural competence, please review the NASW resource: [Standards and Indicators for Cultural Competence in Social Work Practice](#)
5. Need Help? Select the “Help” option in the member dashboard to access FAQs or contact us.

How to print Handouts

On a Mac

- Open PDF in Preview
- Click File
- Click Print
- Click dropdown menu on the right “preview”
- Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.

On a PC

- Open PDF
- Click Print
- Choose # of pages per sheet from dropdown menu
- Choose Black and White from “Color” dropdown

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

How to Navigate Naturally in a Noisy World

Tammy Stender, AuD, CCC-A

Learning Outcomes

- 1) Explain why it is so important for people to navigate as naturally as possible in the diverse and noisy listening situations we encounter on a daily basis.
- 2) Describe the three main technological advancements of ReSound directionality introduced in ReSound OMNIA.
- 3) Detail how a multiband directionality approach provides a myriad of benefits for a user to hear as naturally as possible, for their own voice as well as for all the sounds around them.



One of the biggest innovation challenges the industry faces is:

hearing in noise

BIHIMA Surveys Audiologist Panel on Issues Facing Hearing Industry, Jun 2019



hearing in noise

80%

Of hearing aid
users still
challenged by
hearing in noise*

Experienced users **struggle**
to find solutions that allow
them to participate in
conversations without
compromise or workarounds**

- MarkeTrak 10, 2019
- ** Fidelity Study – Barriers and Triggers, 2021



What Do We Know about Hearing in Noise from Hearing Aid Users?

- Hearing aid performance including hearing in noise was the number 1 driver of hearing aid satisfaction and the biggest unmet need was sound issues including background noise.

(MarkeTrak 10, Picou, 2020)



What Do We Know about Hearing in Noise from Hearing Aid Users?

- Over 2000 current hearing aid users ranked hearing in noise as the most important attribute when purchasing a hearing aid.

(Manchaiah et al., 2021)



What Do We Know about Hearing in Noise from Hearing Aid Users?

- Hearing-in-noise came in last when respondents were asked to rate 17 listening situations with regard to hearing aid satisfaction (n= 536)

(EuroTrack, 2018)



Based on These
Studies, Our Focus is
Naturally...

**Developing Better
Technology for
Hearing-in-Noise**

**To Help Individuals
Navigate Naturally
in a Noisy World**

Introducing ReSound OMNIA



Feels natural

Adapts organically to the person for best wearing comfort all day long and ease of use

Sounds natural

Emulating both the natural hearing process and the natural way we listen as closely as possible, helping the brain to hear as nature intended

Connects users naturally

Inspired by the way we intuitively interact with other devices

Organic Hearing

Sounds natural, feels natural, and connects you naturally to the world



ReSound OMNIA
sounds natural...
But What IS Natural?





Spatial Cue Preservation



Binaural Listening



ReSound GN

Speech Intelligibility

Spatial Cue Preservation Strategy: When You Need to Hear Around You

ReSound GN



Binaural Listening Strategy: When You Want to Focus AND Maintain Awareness

ReSound GN

Speech Intelligibility Strategy: When Understanding is Most Important

ReSound GN



Spatial Cue Preservation



Binaural Listening



Speech Intelligibility

Purpose of ReSound 360 All-Around: Support Natural Auditory Navigational Abilities



360 All-Around



How Does ReSound 360 All-Around Achieve This Goal?



Widening our
sound
environment
mapping




Speeding up our
processing and
steering for faster
mode detection and
selection




**Increasing the
precision** of our
beamforming for our
best speech clarity in
any environment



All working **at the same time across any sound environment**,
enabling users to automatically tune into what they want to listen to
and not have the hearing aids choose for them



**Widening our
sound environment
mapping**




Widening our sound environment mapping

360 All-Around




360 All-Around, when in Binaural Listening mode, streams the off-axis sound from the directional side to the omnidirectional side, providing more auditory awareness in difficult environments (like a popular café)



When Do You Want Multiple Strategies Available?




Café or other complex listening situation




Widening our sound environment mapping





- Compared ReSound ONE All Access Directionality to ReSound OMNIA 360 All-Around
- Audiogram+ targets
- Power domes
- Listening conditions for Binaural Listening mode
- Target: Female Speech at 120 degrees
- Dantale II speech recognition in noise test
 - SRT test to determine 50% correct performance, resulting in a dB SNR score
 - Lower score = better!





Widening our sound environment mapping






Speeding up our
processing and
steering for faster
mode detection
and selection

Speeding up our processing and steering for faster mode detection and selection




New steering system that switches 2X as fast, for clearer distinctions among listening modes




**Increasing the
precision** of our
beamforming for our
best speech clarity in
any environment

New ReSound Front Focus Program –


Selected by
the **Individual**
When They
Want to Be In
the Driver's
Seat!



ReSound Binaural Beamforming Directionality – with Weighting to Help Users Hear Better in Difficult Noisy Situations



Applied in Binaural
Listening Strategy



Applied in Speech
Intelligibility Strategy
and Front Focus
Program






Improvements to the binaural beamformer, providing even more precision for our best speech clarity



- Compared ReSound Ultra Focus (in ReSound ONE) to ReSound Front Focus (in ReSound OMNIA)
- Audiogram+ targets
- Power domes
- Set to ReSound Front Focus or Ultra Focus program
- Target: Speech at 0 degrees
- Dantale II speech recognition in noise test
 - SRT test to determine 50% correct performance, resulting in a dB SNR score
 - Lower score = better!




Although there was little separation of the speech and noise, making it a difficult test, results revealed ceiling effects; both beamformers performed well!

To avoid ceiling effects and to truly validate the beamforming directionality, task was made **more** difficult, by applying a bandpass to the speech and noise from 500-4000 Hz.






Improvements to the binaural beamformer, providing even more precision for our best speech clarity



Omnidirectional

Ultra Focus

Front Focus




+4.25dB
better than omni

+8.61dB
better than omni



Improvements to the binaural beamformer, providing even more precision for our best speech clarity



What is a Doubling of Intensity Equal in Decibels?

60 dB




+

60 dB




= 63 dB



Thus, 4.5 dB represents ~150% better performance.





Improvements to the binaural beamformer, providing even more precision for our best speech clarity, compared to omnidirectional



* Jespersen et al (2022)



Improvements to the binaural beamformer, providing even more precision for our best speech clarity



Better optimized for the head-related transfer function = better precision for better clarity



Ultra Focus

+4.25dB

better than omni



Front Focus

+8.61dB

better than omni



What Is the Difference in Directional Response between 360 All-Around and Front Focus in This Situation?

360 All-Around applies this response according to the environment (i.e., in Speech Intelligibility mode)

Front Focus is user-selected, when the user wants the best SNR!

The ONLY difference may be the Directional Mix setting (crossover frequency between omnidirectional and directional processing)



Navigating Naturally in a Noisy World

ReSound GN

Scorecard

Hearing Others' Speech Well in Noise




Hearing *Your Own* Speech in Noise





OH NO!

Sound Like Yourself



ReSound Multiband Directional Processing



Navigating Naturally in a Noisy World

ReSound GN

Scorecard

Hearing Others' Speech Well in Noise





Hearing *Your Own* Speech in Noise



Comfortable Navigation in Noise



Other Benefits of ReSound Multiband Directional Processing: Sound Quality and Spatial Hearing




Head-Related Directional Characteristics: Impact on Localization Cues

Low frequencies travel around the head with limited intensity changes, in an omnidirectional pattern

- Interaural Timing Difference (ITD) cues, as sound waves travel around the head

High frequencies are affected by head shadow and are more like a directional response

- Interaural Level Difference (ILD) cues, as sound waves are attenuated by head shadow



Wind Noise Mechanics


Hearing aid wind noise is caused by turbulent air flow around the microphone(s)

Results in movement of the mic diaphragm that will be amplified

- Turbulence contains pressure fluctuations, which are converted to electrical fluctuations by the mics

This amplification can overload the hearing aid, especially for **low frequencies**


- Spectral peak of wind noise is around 100 Hz at high wind speeds (Larsson & Olsson, 2004) and tapers off around 3000 Hz



Omnidirectional Processing Helps Control Wind Noise

(Thompson, 2000; Kates, 2008)


- Wind noise is spectrally low-frequency and spatially uncorrelated
- When wind is present, distinct air vortexes are created at each microphone of the hearing instrument
- Each of these vortexes is unique, so the signal at each microphone is unique and thereby uncorrelated
- Combining uncorrelated inputs at the 2 microphones of a directional hearing instrument results in an increase of the signal level



Thompson & Dillon, 2002

- Thus, in a fullband directional hearing instrument (with low-frequency equalization), wind noise is amplified further

The ReSound Multiband Directional Processor



The green arrows indicate that the directional signal of the dual microphones of both hearing aids are weighted according to the hearing situation and added to the output signal in order to generate a stronger directional signal.

Navigating Naturally in a Noisy World

ReSound GN

Scorecard

Hearing Others' Speech Well in Noise



Hearing *Your Own* Speech in Noise



Comfortable Navigation in Noise



ReSound OMNIA

Ready to Play...

Guess
The
Answer!

Write it in the chat!



Which listening mode of 360 All-Around emphasizes naturalness and overall sound quality?



● Spatial Cue Preservation

Which listening mode of 360 All-Around uses the brain's binaural processing to shift focus without losing awareness?




■ Binaural Listening

Which listening mode of 360 All-Around provides the best SNR in noise when speech is in front?



▲ Speech Intelligibility

The proximity effect is...



▲ A longstanding issue with directional microphones

◆ Made even worse when "bass boost" is applied

● Directional mics amplify LF sounds that are near more than far away

What new feature of ReSound OMNIA replaces Ultra Focus?



● Front Focus

What is an advantage of a multiband directional approach?



▲ Better access to localization cues than fullband directionality

◆ Better sound quality as compared to fullband directionality

● More natural own voice perception as compared to fullband directionality

ReSound M&RIE: Also a good option for ReSound OMNIA!



Reduced listening effort

Crystal clear phone calls
& natural phone placement

A more natural hearing experience
through preferred sound quality

Improved localization and depth perception
that improves over time

Increased wind noise reduction
and better sound quality in wind

Ready to hear the difference?

ReSound GN

Q&A

Tammara Stender

tstender@gnresound.com

