

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
3. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.
4. Social Workers: For additional information regarding standards and indicators for cultural competence, please review the NASW resource: [Standards and Indicators for Cultural Competence in Social Work Practice](#)
5. Need Help? Select the “Help” option in the member dashboard to access FAQs or contact us.

How to print Handouts

On a Mac

- Open PDF in Preview
- Click File
- Click Print
- Click dropdown menu on the right “preview”
- Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.

On a PC

- Open PDF
- Click Print
- Choose # of pages per sheet from dropdown menu
- Choose Black and White from “Color” dropdown

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

A solid red square is located in the bottom left corner of the page.

Academia

CLAYTON FISHER

M.Cl.Sc., Audiologist

(Hands On)
**REM
Basics
with CARL**



Wednesday, October 26, 2022

12:00PM EST



Another how-to REM workshop?



Yep, but here's how this one's different:

- 1) Focuses entirely on **practical** aspects of REM, **not on the theoretical**
- 2) We are going to **do a fitting** together, using real equipment, **not just talk** about it
- 3) We are going to **make mistakes** together and try to **troubleshoot** problems

REM

=

Rreal

Ear

Measurement

About Me



I have worked with various organizations to help clinicians do their job better

- Optimizing equipment to make clinician's lives easier, more streamlined
- Creating best-practice protocols that optimize patient outcomes but are not overly burdensome for clinicians
- Training and supporting clinicians to ensure they have the required skills and understanding to perform optimally clinically

Workshop Goals



Learning Objectives:

- 1) Feel more comfortable with your REM **equipment and software**
- 2) Feel more comfortable with **calibration** and **techniques** needed to ensure accuracy
- 3) Feel comfortable **interpreting** your REM results and fitting data and **improving** upon it

...and get motivated to use your equipment!

Workshop Outline



- **Part 1: Context on today's talk (the “Why”)**
 - What we will not be talking about
 - Our relationship with our audiology equipment
 - A little about my *non*-audiological background
- **Part 2: Getting comfortable by *doing***
 - Practical REM equipment and patient set-up
 - Ensuring measurement accuracy during REM
 - Interpreting REM results and making improvements

What we won't cover today




We will not discuss:

- The well-documented, known benefits of REM
- Why choosing to ignore best-practices is unethical

For more info on these topics...

- Check out the videos on my YouTube Channel:



<https://www.youtube.com/channel/UCvaPjC3au95PAUtD6WRI0Hw>

What we won't cover today



I promise that we will not discuss:

- Terminology
- Prescriptive methods

For more info on these topics...

Check out:

- A textbook
- 20Q articles (Gus Mueller)
- AO REM videos

REUR vs REUG

REAR vs REIG

RECD vs REDD

DSL vs NAL

My assumptions



a c a d e m i a

If you're participating today, I assume:

- You know REM is key to ensuring ideal patient outcomes, and
- You want to be doing REM routinely, but you're not, or
- You're doing REM routinely, but you want to understand it better and feel more confident



Why I use REM



academia

“Audibility always wins...” – Gus Mueller

- If you want to ensure audibility is adequate, you need to measure it; improve it
- The *improvement* in audibility achieved in a fitting is a key factor in predicting the degree of benefit with hearing aids
- Thank goodness for automatic aided SII, unaided SII, and SII *difference* calculations

SII (%)

R		L	
A	U	A	U
39	41	63	64
-2		-1	
54	40	67	64
+14		+3	
62	40	70	65
+22		+5	

SII = Speech Intelligibility Index = % of speech audible to patient


Context for this workshop



academia

- **REM systems have evolved**
 - Computer-based
 - User-friendly
 - Portable
- **So has the way we learn, teach and practice**
 - Patient simulators: meet **CARL!**
 - Allows practice of key techniques
 - Allows safe environment for mistakes

The Good:



Software suite

Maestro



AHEAD
SIMULATIONS

Context for this workshop



academia

The *Not-So-Good*:

- **I've realized, through experience, that many clinicians**
 - Are not gear-savvy
 - Are not “technical” per se
 - Are not fully comfortable with their equipment
- **In addition...**
 - Training varies greatly across programs
 - Practicums can help create good habits, or bad ones
 - It is easiest to just maintain the status quo

Myth Buster



a c a d e m i a

Which of these do you think is the biggest barrier to REM uptake?

- a) Ignorance of benefits
- b) Cost of equipment
- c) Access to equipment
- d) None of the above



Myth Buster




academia

According to a 2010 survey of clinicians by The Hearing Journal, about half (48%) of clinicians who have a REM system in their office

don't use it!

Anecdotaly, I can confirm this is true from working for two larger hearing aid clinics



Why then?

Mueller HG, Picou EM. Survey examines popularity of real-ear probe-microphone measures. Hear Jour. 2010;63(5):27-32.



The Big Barrier



academia

Lack of training and familiarity with equipment and procedures greatly impacts the use of REM and thus patient outcomes

Clinicians who don't use REM are often not comfortable with:

- Hardware and software
- Calibration and ensuring accuracy
- Interpreting results and making adjustments to improve the fitting

Let's “Talk Shop”



academia

“A piece of gear you can depend on is honestly a big part of my mental and emotional security...”

- Brian Lam, Wirecutter



Our Relationship with Gear



academia

Obviously, we're not *all* gear-heads...

Where do you place yourself on this scale?



**I only have / want the minimum
knowledge I need to get by**

**I want to know ever feature,
usage and specification**

Note: I'm probably around here

A red arrow originates from the text 'Note: I'm probably around here' and points to the scale bar, specifically towards the number 5.

My background outside Audiology



a c a d e m i a

I was an aspiring musician before becoming an audiologist

- yep, I played guitar in heavy rock band

As such, I've always dealt with equipment

- Guitars, amps, pedals, patch cords, etc.
- Troubleshooting on stage: feedback, intermittent sound, no sound, etc.



Lean In... to your equipment



Equipment is part of audiology... so

- Ask a more experienced colleague to give you an overview of equipment you're unfamiliar with
- Use a cancelled appointment as an opportunity to learn a new feature, or troubleshoot issues that came up in clinic earlier that day or week
- Take a little time over lunch to fool-around with your equipment, look over the user guide, etc.
- Be curious and have fun... Your REM system should be your best friend! Get to know it.

...and be thankful you don't have to use legacy equipment, as pictured below



Disclosure



a c a d e m i a

Truthfully, I have not been using Inventis equipment for that long...

- I was using a different manufacturer, and have switched
 - Each brand of device is different
 - Change is hard at first!
 - I'm still learning... but the best way to learn is by *doing*!




Part 2: Let's fit CARL!



CARL is a patient simulator with **realistic ear canal acoustics**

- No patient needed, no risk of harm
- He's great for practicing common clinical procedures:
 - REM
 - Wax removal
 - Ear impressions
 - Etc.




<https://aheadsolutions.com/>



CARL's Audiogram



academia



	R	L	BIL
AC	○	×	⊗
AC Mask.	△	□	△
UCL	u	u	
BC	∩	∩	
BC Mask.	∩	∩	
FF	≡	≡	≡
FF Mask.	≡	≡	≡
FF Aided	A	A	A

Audibility index		
	R	L
AC	64 %	42 %
Pure tones average		
	R	L
AC	25	30
BC	-	-
Freq. [Hz]: 500, 1k, 2k		



<https://aheadsims.com/>

Step 1: Patient Set-up



academia

After positioning patient, perform

1) Speaker (Room) Equalization

- Modifies the speaker frequency response, based on acoustics of room
- Weekly, daily, or anytime you change something or are in a new space

2) Probe Tube Calibration

- Makes resonances of probe tubes acoustically invisible
- Also, probably fine to do once a week

A screenshot of the 'REM PREFERENCES' settings screen. The screen has a white background with a light gray sidebar on the left containing two menu items: 'MEASUREMENTS SETTINGS' and 'NOTIFICATIONS SETTINGS'. The 'NOTIFICATIONS SETTINGS' menu item is highlighted with an orange border. To the right of the sidebar, there are several settings. At the top right, there are three buttons: 'NEW PATIENT', 'DAY', and 'WEEK'. Below these buttons, there are five toggle switches. The first three are currently turned off (gray), and the last two are turned on (orange).

REM PREFERENCES	
MEASUREMENTS SETTINGS	Equalization expires every
NOTIFICATIONS SETTINGS	Enable speaker equalization reminder
	Enable probe calibration reminder
	Check equalization at startup
	Enable ambient noise warning
	Check differences between sides

Patient Set-up



academia

3) Probe tube Insertion

- Use lanyard to hold tube in place
- Depth matters; go deep, or go home
- Don't be shy, you can't perforate an eardrum with a probe tube!
- Check probe placement with otoscopy
- Luckily, CARL does not have wax, bendy canals, hairy canals, thin canals, etc.

After calibrating...



Step 2: Software Set-up



- **Check your fitting parameters**
 - Choose your prescriptive targets
 - HA style
 - acoustics, etc.

DSL-V5 SETTINGS

PATIENT DATA		
AUDIOLOGICAL INPUT		
HEARING AID		

Number of hearing aids	Binaural ▼
Type of hearing aid	BTE ▼
Number of channels	16 ▼
Circuit type	WDRC ▼
Venting	Open ▼

DSL-V5 SETTINGS

PATIENT DATA		
AUDIOLOGICAL INPUT		
HEARING AID		

Transducer	Insert	TIP MOLD
ABR	☒	NHL EHL
Use BC		YES NO
Applied UCL		PREDICTED MEASURED
RECD		PREDICTED MEASURED
REUG		PREDICTED MEASURED
RECD type		HA1 tip ▼
REUG type		0° ▼

SAVE CLOSE

Optional: Run Resonance Curves



academia

1) Check ear canal resonance with pink noise

- REUR should look reasonable

2) Check occluded response

- Mute HA, insert, and run REOR
- How much does HA act like an earplug when not amplifying?

3) Troubleshoot, if needed

- Has probe tube moved?



Step 3: Aided Speech



academia

Next, “best-fit” your hearing aids using your favorite fitting formula: DSL, NAL, Manufacturer

- **Complete your level calibration**
 - If the fitting is open, or heavily vented, choose “open fit” and mute HA for cal.
- **Start your aided speech measurements**
 - Despite unaided SII predictions, it is still a good idea to run an unaided speech curve at 65dB



Step 4: Interpret Results



academia

Since **we've dotted our i's** and **crossed our t's**, we are confident that the picture that is painted is a **true reflection of audibility in the patient's ear**

- What does the fit-to target look like?
- How much improvement in audibility is there?
- Can you improve audibility while keeping average speech "Comfortable, but Slightly Loud"



+12

+23


Step 5: Adjustments



academia

Make adjustments to target to improve audibility, **prioritizing 2 - 4 kHz** for open fittings and **0.5 – 4 kHz** for closed fittings

- **Use “On-top” mode**
 - Facilitates adjustments
- **Again, complete your level calibration**
 - If the fitting is open, or heavily vented, choose “open fit” and mute HAs for cal.



Step 6: Run Other Inputs



academia

Deviation from target should be examined for different input levels

- **Soft speech** (55 dB) should be “Comfortable, but **Slightly Soft**”
- **Loud speech** (75 dB) should be “**Loud, but OK**”
- **MPO** should be kept **below UCLs** especially in ear canal resonance for open fittings, and
 - Not too low for closed fittings



Troubleshooting



academia

What common problems issues can you run into?

Probe tube issues:

- Shallow probe tube placement **depth**
 - Can move out, especially if HA is removed
- Probe tube gets blocked with **wax**
 - It only takes a tiny amount of soft wax to block a probe tube completely



Making Mistakes



academia

We all make clinical mistakes; the goal is to only make them once.

- **Errors with REM**
 - Swapping reference mics
- **Errors with HAs**
 - Wrong receiver strength in software
 - Swapped receiver sides






Catching Mistakes



academia

So, we all make mistakes...
wouldn't it be great if we
could catch them **while**
they're happening?

- REM shows you problems that you would not otherwise know about



Summary



a c a d e m i a

- Real-Ear Measurements are essential for improving patient outcomes.
- REM should not be viewed as being overly technically difficult; it's not that hard!
- To feel more comfortable, we need lean into our gear, get to know it, roll up our sleeves, and try things out.
- CARL is great a patient to have in the office to practice on.
- Start a friendship with your REM system; it will last a lifetime!



More about CARL



a c a d e m i a

- Today we were using “Base CARL”
- There is also “**Camera CARL**” which gives you visual feedback about where you are in the ear canal
- There is also **Audio CARL** who responds to acoustic stimuli...



Thank you for your attention.

— QUESTIONS?

Contact: info@treatyourhearing.ca

www.inventis.it

