

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
3. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.
4. Social Workers: For additional information regarding standards and indicators for cultural competence, please review the NASW resource: [Standards and Indicators for Cultural Competence in Social Work Practice](#)
5. Need Help? Select the “Help” option in the member dashboard to access FAQs or contact us.

How to print Handouts

On a Mac

- Open PDF in Preview
- Click File
- Click Print
- Click dropdown menu on the right “preview”
- Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.

On a PC

- Open PDF
- Click Print
- Choose # of pages per sheet from dropdown menu
- Choose Black and White from “Color” dropdown

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

A solid red square is located in the bottom left corner of the page.

Real-Ear Measures

A Deep Dive

AMANDA WOLFE, Au.D. EDUCATION AND TRAINING AUDIOLOGIST



Hear better. Live better.

Learning Outcomes

1

Participants will be able to explain the benefits of utilizing real-ear measurements in their practice

2

Describe how to prepare for real-ear measurements

3

Describe the real-ear measurement process



Hear better. Live better.

Best Practice



Contributes to higher overall satisfaction with hearing aids and their hearing care professional, as well as an improved patient experience¹.



Includes individual ear information for a **more precise, personalized hearing aid fitting**.



The best way to verify the performance of a hearing aid² according to best practices.

1. Amlani AM, Pumford J, Gessling E. Improving patient perception of clinical services through real-ear measurements. *Hearing Review*. 2016;23(12):12-21.

2. Valente M, Abrams H, Benson D, et al. Guidelines for the Audiological Management of Adult Hearing Impairment. *Audiology Today*. 2006;Vol 18.

Goals of REM



Hear better. Live better.

Goals



Verify
Audibility



Verify
Output



Ensure Patient
Benefit



Troubleshoot
Complaints

Soft speech and
environmental sounds
are audible



Average speech
is comfortable



Loud sounds
are not uncomfortable



Independent Fitting formulas

NAL-NL1 (*Byrne et al, 2001*)

- Maximizes speech intelligibility
- Keeps overall loudness at a level no greater than that perceived by a normal-hearing person

NAL-NL2 (*Keidser et al, 2011*)

- Revised version of NAL-NL1, but provides a different gain-frequency response shape
- Slightly higher compression ratios for mild/moderate hearing loss

DSL v5 (*Beck et al, 2007*)

- Provide maximum audibility
- Maintain comfortable loudness

NAL-R (*Byrne & Dillon, 1986*)

- Intended for linear hearing aids
- Amplifies speech so that all frequency bands contribute equally to its loudness



Integrating REM into your Clinical Practice



Hear better. Live better.



Hearing Test



Fitting



Follow-up



Hearing Test



Fitting




Follow-up



REM Terminology



Hear better. Live better.

An underwater photograph showing sunlight rays filtering down from the surface. The water is a deep blue-green. At the bottom, there is a dark, textured coral reef. Several white text labels are overlaid on the image.

REUR

REOR

REAR

REOG

REIG

REAG

RECD

SII

Response

An *absolute* measure of SPL

$$\text{Input} + \text{Gain} = \text{Response/Output}$$

Gain

A response that is expressed *relative* to a specified input or reference response

$$\text{Output} - \text{Input} = \text{Gain}$$

REAR Real-Ear Aided Response

(ANSI S3.46-1997)

The SPL, as a function of frequency, at a specified measurement point in the ear canal for a specified sound field, with the hearing aid (and its acoustic coupling) in place and turned on.

The frequency response of a hearing aid that is turned on, measured in the ear canal, for a particular input.

REAG Real-Ear Aided Gain

(ANSI S3.46-1997)

Difference in decibels, as a function of frequency, between the SPL at a specified measurement point in the ear canal and the SPL at the field reference point, for a specified sound field, with the hearing aid (and its acoustic coupling) in place and turned on.

The gain of the hearing instrument across frequencies, measured in the ear canal. Subtract the input signal level from the REAR across frequencies to obtain the REAG.

$$\text{REAR} - \text{INPUT} = \text{REAG}$$

SII Speech Intelligibility Index

The SII score is a way of quantifying the likely intelligibility of speech based on the proportion of speech cues that are audible to the listener.

It is a measure ranging between 0.0 and 1.0 “that is highly correlated with the intelligibility of speech.” (ANSI, S3.5, 1997, p. 1).



On-Ear Measures Setup

Audioscan Verifit 2



Hear better. Live better.

Audioscan Verifit

Speechmap





Set-up Steps

- ① Calibration
- ② Enter audiogram
- ③ Select fitting formula
- ④ Select device type/coupling
- ⑤ Select stimulus
- ⑥ Patient placement
- ⑦ On-ear probe tube placement
- ⑧ Hearing aid set-up
- ⑨ Patient preparation

SETUP

Calibration


Compares frequency response of probe tube to flat frequency response of the external reference mic, and compensates for this difference in all subsequent measurements



STEP 1

Calibration

Fold the probe tube over and place the probe tube opening on top of the reference microphone.



STEP 1

Calibration

To Calibrate

- Hold the probe dock 6-36 inches from the speaker
- Select Calibrate




STEP 1

Calibration

To Calibrate

- View calibration curve
- Repeat for other probe dock




Enter Audiogram

- Audiometric entry**
- | 250 | 500 | 750 | 1k | 1k5 | 2k | 3k | 4k | 6k | 8k | 10k | 12k5 | | | |
|-----|-----|-----|----|-----|----|----|----|----|----|-----|------|---|-----------|---|
| 50 | 50 | | 50 | | 55 | | 60 | | 65 | | | X | Threshold | X |
| | | | | | | | | | | | | X | BCT | X |
| | | | | | | | | | | | | X | UCL | X |
- 4d 15
-
- L**
- | 250 | 500 | 750 | 1k | 1k5 | 2k | 3k | 4k | 6k | 8k | 10k | 12k5 |
|-----|-----|-----|----|-----|----|----|----|----|----|-----|------|
| 50 | 50 | | 50 | | 55 | | 60 | | 65 | | |
- R**
- | 250 | 500 | 750 | 1k | 1k5 | 2k | 3k | 4k | 6k | 8k | 10k | 12k5 |
|-----|-----|-----|----|-----|----|----|----|----|----|-----|------|
| 50 | 50 | | 50 | | 55 | | 60 | | 65 | | |

STEP 3

Select Fitting Formula



- Determines fitting targets based on fitting formula philosophy
- Match fitting formula selected in the fitting software
- Approximate fitting formula based on manufacturer proprietary formula*



STEP 4

Select Device Type and Coupling/Vent

- Ensures appropriate target generation
- Determines need for sound field equalization



Sound Field Equalization

Concurrent

- For occluded fittings
- Reference microphone is on during measurements
- If the patient moves during the measurement, the system will adjust the signal level as needed.
- May result in an underestimation of hearing aid gain if used with open fit.

Stored

- For open fittings
- Reference microphone is off during measurements
- Patient must remain still during the measurement or equalization should be repeated if patient moves.
- Will prevent underestimation of hearing aid gain due to leakage.

STEP 5


Select Stimulus

ISTS

- Female talker
- Segments of 6 languages sliced together to match the average speech spectrum

Speech-std

- Carrot passage
- Filtered to provide the long-term average speech spectrum (LTASS)




STEP 6

Patient Placement

To Position

- Choose a quiet location
- Position patient 18-24 inches in front of the speaker




*Follow instructions from real-ear system manufacturer

STEP 7

On-Ear Probe Tube Placement

To Position

- Clip probe dock to patient's shirt
- Hang ear loops over pinna and adjust to fit near the head and away from acoustic influence of pinna or hearing aid



STEP 7

On-Ear Probe Tube Placement

To Position

- Positioning the probe tube
 - **Goal:** Within 2-5mm of eardrum
 - Adult males: ~30mm
 - Adult females: ~28mm
 - Children: ~20-25mm



CAUTION: Do not poke the eardrum with the probe tube

STEP 8

Hearing Aid Setup



STEP 9

Patient Preparation

- Remain quiet and still
- Will be listening to speech at different volume levels
- Do not need to repeat back or even understand stimulus
- If signal gets uncomfortable loud, raise hand





Hearing Test



Fitting



Follow-up



Verifit 2

Advanced Verification

On-ear



Speechmap

Directional

Noise reduction

Manual control

Calibration

Test box



WRECD

Occlusion

Feedback

Insertion gain

ANSI / IEC

Distortion

Input/Output

Battery drain

Multicurve

Telecoil



Left/Right



Print



Session



Help



Setup

Right-click mouse at any time
to return to this menu.


Version: 4.22.2

August 7, 2020

Serial number: C1776

Network address: 10.33.16.123






- ① Audiometry Entry
- ② Start test
- ③ Aided envelope
- ④ Aided SII
- ⑤ Unaided predicted SII


Max TM SPL 125

L



Max TM SPL 125

R




250 500 1k 2k 4k 8k 16k


On-ear Graph Dual view SII L R Adult Insert + foam Non tonal SPL

RITE Occluding Audiometry RECD Average BCT N/A Bin Yes Counseling	Speech-ISTS	0		0	Speech-ISTS
	Soft (50)	SII	1	SII	Soft (50)
	Speech-ISTS	10		10	Speech-ISTS
	Avg (65)	SII	2	SII	Avg (65)
	Speech-ISTS	56		56	Speech-ISTS
	Load (80)	SII	3	SII	Load (80)
	MPQ	N/A		N/A	MPQ
	BS	N/A		N/A	BS


RITE
 Occluding
 Audiometry
 RECD Average
 BCT N/A
 Bin Yes
 Counseling




- ① Audiometry Entry
- ② Start test
- ③ Aided envelope



SII (%)	% Cor. Digits	% Cor. IEEE Sentences	% Cor. NU6 Words
0	0	0	0
5	30	4	4
10	80	9	8
15	96	13	11
20	98	19	19
25	99	43	27
30	100	60	36
35	100	74	45
40	100	89	56
45	100	94	63
50	100	96	70
55	100	98	76
60	100	99	81
70	100	99	87
80	100	100	93
90	100	100	96





① DSL REAR90 Target

② Estimated UCL

③ Measured MPO

④ LTASS 65

⑤ LTASS 50

⑥ SPL Threshold


Max TM SPL 125

L

Max TM SPL 125

R

Pre-Adjustment



Post-Adjustment




Using SII for Counseling

	50 dB SPL		65 dB SPL		80 dB SPL	
	Left	Right	Left	Right	Left	Right
Unaided SII	0 %	0 %	10 %	10 %	56 %	56 %
Aided SII (pre)	6 %	11 %	32 %	34 %	63 %	64 %
Aided SII (post)	31 %	35 %	54 %	57 %	69 %	70 %



Hear better. Live better.



Match to Target (S_{II}) 55

Patient Preference (S_{II}) 22




- ① Encourage acclimatization to higher gain settings and provide a volume control
- ② Adjust hearing aid gain to a middle ground
- ③ Program hearing aid memories with different gain settings
- ④ Start at lower gain settings and increase hearing aid gain over time





Hearing Test



Fitting




Follow-up

Now What?


Audibility and comfort verified

Our Role


Ensure patients are satisfied
with their hearing aids



Hear better. Live better.



Real World Balancing Act



Hear better. Live better.

Greater perceptual
benefit with REM

Reduced Listening Effort

Improved Speech Understanding

Reduced Patient Visits





Real-Ear Measurements

Equipment Troubleshooting



Hear better. Live better.

Minimal to No Aided Benefit

Culprit

Compromised Probe Tube

Hearing aid receiver blocked

Earmold outlet blocked

Troubleshooting Tips



Check the probe tube

- Clear of cerumen
- Positioning/depth
- Replace the probe tube



Check the hearing aid or earmold

- Clear of cerumen
- Listening check
- Test box measurements



Check the REM system and hearing aid setup

- Demographics, audiogram, acoustic coupling
- Fitting formulas match



Hear better. Live better.

Curves are not Changing

Culprit

Disconnected hearing aid

Compromised probe tube

Hearing aid receiver blocked

Earmold outlet blocked

Troubleshooting Tips

- ✓ Check the connection
- ✓ Check the probe tube
 - Positioning/depth
 - Replace the probe tube
- ✓ Confirm adjustments are being made correctly
- ✓ Check the hearing aid or earmold
 - Clear of cerumen
 - Listening check
 - Test box measurements



Hear better. Live better.

Verifying audibility

STARTING POINT



Real-world
Performance

Sound
Quality

Validation
Measures

Learn More

American Academy
of Audiology

<https://www.audiology.org>

American Speech and
Hearing Association

<https://www.asha.org>

Audioscan Verifit[®]
User's Guide 4.28

**Visit audiologyonline.com and [StarkeyPro.com](https://starkeypro.com)
to register for more training classes**

Thank you



Hear better. Live better.