

The Prescriptive Hearing Aid Fitting Formula Landscape

ANDREA HANNAN DAWKES SENIOR AUDIOLOGIST EDUCATION & TRAINING



Hear better. Live better.

The Prescriptive Hearing Aid Fitting Formula Landscape

How to Ask a Question

Click Q & A, type your question, click Send

Where to get the Handout

A PDF handout for this presentation may be accessed through your Pending Courses on Audiology Online or by clicking on the link in Chat within the classroom

How to Earn CEUs

CEUs are available with a paid Audiology Online Unlimited CEU Access membership. After completing the course, go to your Pending Courses on Audiology Online and pass the multiple-choice exam

Assistance or Tech Support



Use the Q&A window



Call Audiology Online | 800-753-2160



customerexperience@continued.com

Learning Outcomes

As a result of this course, participants will be able to:

List at least three different goals of hearing aid fitting formulas.

Identify the differences between generic and proprietary hearing aid fitting rationales.

Describe the features of at least three hearing aid fitting formulas and the patient populations for whom they are intended.



Agenda

History of Fitting Formulas

Independent/Generic Formulas

Proprietary Formulas

Clinical Considerations

“The goal in hearing aid selection and fitting is to provide comfortable hearing instruments that will meet the lifestyle needs of the patient while easing the patient’s communication difficulties and maximizing performance in different listening environments.”

J. L Northern, 2011

Hearing Aid Selection & Fitting

Product Family
Technology Tier
Style
Matrix
Coupling
Venting
Advanced Features

Hearing Test Results
(AC, BC, MCL, UCL)

Preferences

Lifestyle

Unique Listening Needs & Goals

Expectations

Performance in Quiet & Noise

Previous Experiences

Cognition/Memory Screening

Physical Dexterity



Hearing Aid Satisfaction

HA Performance
and sound

0.67

HCP Effectiveness

0.42

HA Physical Qualities

0.25

HA Maintenance

0.24

Cost
(upfront and ongoing)

0.20

MarkeTrak 10 / Hearing Industries Association, 2019

Top Drivers of Hearing Aid Performance and Sound

MarkeTrak 10 / Hearing Industries Association, 2019



Natural Sounding

Rich and Clear



Minimize background noise



Comfort for loud sounds



Ability to tell direction

Goals of Hearing Aid Fitting Rationales



Minimize initial **rejection**



Maximize **speech intelligibility**



Maximize usage of residual **dynamic range**



Maximize **patient comfort**



Maximize **sound quality**



Maximize patient **satisfaction**

POGO

NAL-NL2 Tonal

Manufacturer
Proprietary

Half Gain Rule

DSL v5.0 Pediatric

NAL-RP

FIG6

NAL-R

BERGER

Fitting Formulas

VIOLA

NAL-NL2 Pediatric

DSL [i/o]

Pascoe

DSL v5.0 Adult

LIBBY 1/3

NAL-NL1

NAL-NL2

POGO II

Hearing Aid Fitting

“About 27,000,000”

Hearing Aid Fitting Formula

“About 5,740,000”





1802
Floral
Aurolase
Phone



1808
Ear Horns



1925
Carbon
Amplifiers



1938
Vacuum Tube
Hearing Aids



1944
Lybarger
Half Gain Rule

1944
Transistors




1946
Harvard Report

1946
Body Aids



1955
Early BTE
and Custom
hearing aids



1

9

4

0


Watson and Knudsen

Selective amplification

Frequency response specific to degree of hearing loss and audiometric configuration

Suprathreshold method for determining required gain and frequency response based on the patient's most comfortable loudness (MCL) levels across frequencies

First prescriptive fitting formula



1

9

4

6

Harvard Report


Davis, et al.

Common frequency response

Most patients performed best with the same frequency response

Flat or 6dB/octave rising response between 300 and 4000 Hz

Primary objective of hearing aid evaluations was the “intelligibility of speech” using word tests



1

9

4


6

Comparative Approach

Patient is tested with several hearing aids using a variety of speech recognition tasks

Hearing aids could be differentiated from one another if the word recognition scores differed by 8% or more

Primary method for choosing hearing aids until the early 1980s



• 1975

Pascoe (1975)

Frequency-gain characteristics of hearing aids are important for optimizing speech intelligibility with sensorineural hearing loss

• 1978

Thornton & Raffin (1978)

Northwestern University (NU)-6 word recognition task had very poor **test-retest reliability which could impact the ability to identify differences between hearing aids**

• 1980

Studebaker & Hochberg (1980)

“Thirty years of extensive research have failed to demonstrate that hearing aids can be chosen reliably or validly by [the comparative] method in a reasonable period of time.”

• 1983

Walden, et al. (1983)

A comparative procedure was not valuable unless hearing aids with very dissimilar frequency responses were being tested - unlikely in clinical practice

“Real ear measures came on the scene and audiologists began to appreciate the difference between 2cc coupler output and measurements in the real ear along with the impact of earmolds on the hearing aid response” (Killion, 1981)

Today

DSL v5.0 Adult

NAL-RP

BERGER

NAL-NL2

DSL v5.0 Pediatric

Manufacturer
Proprietary

NAL-R

NAL-NL1

NAL-NL2 Pediatric

Prescriptive Fitting Formulas

1

Theoretical
Rationale

2

Amount and
Type of
Audiometric
Information

3

Type of Signal
Processing



Theoretical Rationales

LOUDNESS EQUALIZATION

Loudness across frequency bands is equalized so all are perceived as equally loud

Patient should **perceive overall speech levels similarly to a person with normal hearing**

LOUDNESS NORMALIZATION

Patient should **judge soft, normal and loud sounds like someone with normal hearing**

Maintain normal relative loudness across frequency bands but with **lower frequencies perceived louder - vowels are more intense than consonants**

Amount and Type of Audiometric Information

AUDIOMETRIC INFORMATION


Hearing thresholds (AC, BC)

Transducer (Insert earphones, TDH, Sound field)

UCL / LDL

MCL

Hearing Loss Configuration



Type of Signal Processing

LINEAR

- Gain fixed across input levels
- Gain varies as a function of frequency
- A single gain target is produced

NON-LINEAR


- Gain varies as a function of frequency and input level
- Developed for use with compression hearing instruments to keep sounds within the user's dynamic range
- Multiple gain targets are produced


Formula Fun Fact

PYTHAGOREAN THEOREM

The square of the hypotenuse of a right triangle is equal to the sum of the squares of the other two sides.

$$a^2 + b^2 = c^2$$






Independent/Generic Fitting Formulas

Linear Fitting Formulas

One frequency/gain target is produced



Half Gain Rule

Berger

POGO

Libby 1/3

NAL-R & NAL RP

DSL v5.0

1
9
4
4

Half Gain Rule

Lybarger

Aided thresholds improved by ***half the degree of hearing loss***


Average ***conversational speech will be amplified to the middle of the listener's dynamic range***

Less than half-gain is provided for the very lowest frequency to ***reduce the upward spread of masking***

Basis for subsequent linear fitting methods



Phonemic Classification



Berger et al.


More than half gain provided at 1000 Hz and 2000 Hz - the most important frequencies re: speech cues

Keep less than half gain for the low frequencies to reduce the upward spread of masking and optimize speech intelligibility

Guidance for selecting maximum output (SSPL90) – listener's UCL minus 3 to 6 dB

10 dB of reserve gain

May be used with linear and non-linear hearing aids



1
9
8
6

NAL-R

Byrne and Dillon

Adjacent *speech frequencies should sound equally loud* at a comfortable listening level

Reduce the intensity of louder speech sounds so *all speech sounds are equally loud* when aided

Less high frequency gain than Berger *for steep, high frequency hearing losses*

Target gain looks 'flatter' than those prescribed by most other linear fitting methods

15 dB reserve gain



NAL-RP

Byrne, Parkinson, & Newall

Considerations for severe-to-profound hearing losses

- Amplify less at frequencies where there is severe-to-profound loss
- Amplify more at frequencies where there is better hearing


Cochlear dead regions



Non-linear Fitting Formulas

Used with hearing aids offering compression
so aided response is within the dynamic range

Gain varies as a function of frequency and
input level - gain for soft input levels greater
than gain for moderate or loud inputs levels



NAL-NL1 & NAL-NL2

DSL v5.0

FIG6

VIOLA/IHAFF

NAL-NL1

National Acoustic Laboratories' Nonlinear Fitting Procedure, Version 1

Compression-based method from the linear NAL-R method

Equalize loudness relationships
among the speech frequencies

Gain for several input levels for maximal speech intelligibility

Overall loudness perception of
aided speech is consistent with
that experienced by those with
normal hearing

Less prescribed gain in regions
of severe hearing loss

Includes ***age-appropriate correction
factors*** (e.g., for RECD)

Available for ***pediatric population***



2

0

0

5

DSL v5.0

Desired Sensation Level (DSL) Method
Scollie, et al.

Hearing aid prescription for *infants, children and adults*

Higher listening level for *pediatric patients* than for *adult patients*

Calculates a *full range of targets for any audiogram of two or more thresholds*

Compression calculated with consideration for the *channel architecture*

Targets for aided speech are limited so that *speech peaks do not exceed the listener's UCL*

Full audibility of speech in quiet / less gain at speech frequencies of lower importance with a raised CT in noise

Binaural correction reduces speech targets by 3 dB

Higher targets for *conductive and mixed hearing loss*

2

0

0

8

NAL-NL2

NAL's second generation prescription for fitting WDRC instruments

Attributes of the hearing aid user are considered: ***gender, experience level, binaural listening, age, and language***

Female hearing aid users prefer lower overall gain than male users

New hearing aid users with more than a mild hearing loss prefer less gain with increasing degree of hearing loss than experienced hearing aid users

Unilaterally and bilaterally fitted hearing aid users prefer overall gain levels that vary less than estimated by the bilateral correction factor



2

0

0

8

NAL-NL2

Adults prefer lower overall gain than children

People with severe/profound hearing loss prefer lower compression ratios than predicted when fitted with fast-acting compression.

Slightly more low-frequency gain is prescribed for **speakers of tonal languages** (e.g. Mandarin)




Formula Fun Fact

FUNDAMENTAL THEOREM OF CALCULUS

Establishes the relationship between differentiation and integration and provides a way to evaluate definite integrals without using Riemann sums or calculating areas.

$$\frac{df}{dt} = \lim_{h \rightarrow 0} \frac{f(t+h) - f(t)}{h}$$



Proprietary Fitting Formulas

Developed by hearing aid manufacturers

Linear and non-Linear prescriptive
approaches


Designed to optimize hearing aid
performance based on the unique
compression architecture and signal
processing available



Manufacturer Proprietary Formulas

- Compression type
- Compression speed
- Compression attack and release times
- Number of channels
- Microphone location effects
- Acoustic effects of hearing aid style and coupling





Fitting
Formulas
Generic Non-Linear

Music


Hearing loss impacts the enjoyment of music

Individual preferences for music exist

Music and speech are very different signals



Music vs Speech




Music

Evidenced-based prescription optimized for music

Restores audibility for soft music

Restores desired loudness for loud music

Unique compression TKs & CRs




Dual Speed Compression

Uses both fast-acting and slow-acting compression

Two independent compressors active at the same time

Seamlessly adapts to dynamic situations




Compensating for Loss of Function

Outer Hair Cells

Fast Compression

Inner Hair Cells

Slow compression



2
0
2
3

e-STAT 2.0 Proprietary Fitting Formula

Designed to complement the
Genesis AI compression architecture

*Varies with the audiogram and acoustic stimuli

**Varies with audiogram and with frequency

~4-6 dB more gain
for soft sounds*

10% higher
speech intelligibility scores
for soft speech vs. Evolv AI

More gain at
higher frequencies**
(> 3kHz)

Improved vent modeling

Clinical Considerations


Patient Populations

First Fit Strategy

Impact of Acoustic Options

Advanced Features

Sound Quality



Sensorineural Hearing Loss

Reduced dynamic range

Compression is usually necessary to fit
signals into the patient's dynamic range

Non-linear fitting formula usually optimal

Proprietary / NAL-NL2 / DSL v5.0



Conductive/ Mixed Hearing Loss

Usually want more gain than someone with SNHL

No reduction in dynamic range

Linear processing may be more appropriate

Berger / NAL-R / DSL v5.0



Children

Size

Ear canal properties

Middle ear properties

Central auditory system


Test equipment / procedures

Audibility of speech is crucial – learning
speech and language

DSL v5.0 or NAL-NL2




First Fit Strategy




Audibility / Comfort

Impacts the frequency response

New User



Experienced User



Real Ear Measurements



Best practice by the American Speech-Language-Hearing Association (ASHA) and the American Academy of Audiology (AAA)

Consumer Reports advocates for this practice in their “Hearing Aid Buying Guide”

Predicted real ear responses should not be considered a replacement for the real thing



Impact of Acoustic Options



CLASSICAL MODIFICATION AND EFFECT CHART

	Modification	Effect on Low Frequencies (Below 750 Hz)	Effect on Frequencies Between 750 and 1500 Hz	Effect on Frequencies Between 1500 and 3000 HZ	Effect on High Frequencies (Above 3000 Hz)
TUBING DIAMETER	Larger I.D. Tubing & Horn Tubing	Negligible	Moves peak to higher frequency	Increases height of peak and moves to higher frequency	Increases
	Smaller I.D. Tubing	May reduce below 1 KHz	Moves peak to lower frequency	Reduces height of peak and moves to lower frequency	Large reduction
TUBING LENGTH	Longer Tubing	Increases	Moves peak to lower frequency	Moves peak to lower frequency	Negligible
	Shorter Tubing	Slightly decreases	Moves peak to higher frequency	Moves peak to higher frequency	Negligible



Advanced Features

Feedback management
and initialization

Frequency Lowering

On Demand Adaptive Tuning

*May impact the frequency response
of the amplified signal*

Ski Slope Hearing Loss

More complex than just restoring audibility –
'hearing' the high frequencies may not be useful

Focus is on the mid-frequency region that patients
are relying on for speech understanding

Frequency lowering may be helpful for speech
intelligibility





Sound Quality

Large individual differences

Frequency-gain responses for restoring audibility may result in unacceptable sound quality perceptions

High frequency amplification may be detrimental to speech intelligibility for some users

Sound quality and distortion must be considered

Fitting formulas provide a foundation

Response adjustments allow us to address the unique needs and preferences of each listener

“Individualized adjustments may be needed to obtain a good fit. In this respect, the use of....prescriptive formulas....should probably be regarded as no more than a first estimate of the best frequency-gain characteristic for an individual patient.”

Bratt & Sammeth, 1991

Validation Measures

Optimal performance and sound quality may require adjusting the response away from the prescriptive target

Noise reduction will impact gain at one or more target levels

Various subjective and objective outcome measures are available



Fitting Formula is the Foundation of the Hearing Aid Fitting

Consider type, degree, and configuration of hearing loss, experience, age, and other special considerations

Successful hearing aid fittings strike a balance among audibility, loudness, and sound quality

Both proprietary and generic fitting formulas are valuable

Real ear and outcome measures are important



Thank You.



Hear better. Live better.