

Beltone Solus Max™ 2.0


Developed for Your Expertise and Care

Fitting Flow




Great First Fit Experience

A great fitting experience is a combination of a good foundation and customization. **Beltone Solus Max** provides you the ability to set the basis of your fitting by changing default settings with preferences. Preselecting defaults and experience levels can save time within a fitting session by beginning a fitting with preferred settings. Once the base fitting is done you can make additional changes to tailor the fitting to your patient's unique needs.



Preset settings in Beltone Solus Max

You can access the preference menu of the Beltone Solus Max fitting software in the “Preferences” tab of the home screen.



1 Set first fit sessions launch to Patient Profile Screen

Reason: Starting on this page allows customizations to the initial fit based on previous hearing aid experience and tolerance for amplification.

2 Vent Correction view by default

Reason: The Vent Correction view is a simulation of the vent effect of the selected earpiece from the Physical Properties screen. Selecting ON will result in the display of the low frequency loss with an open or semi open earpiece. This can help determine the appropriate earpiece needed to meet target. Selecting OFF will display the hearing aid amplification regardless of earpiece.

3 Remember to **Save** changes before exiting the tab. Changes will take place upon re-launching the software.



Initial Fit Set Up with Prescribed First Fit

The initial settings for the gain prescriptions appear in the upper left corner of the Patient Profile screen. Select the desired fitting rule and customize by setting the Experience Level and Gain Level % based on the patient's previous hearing aid experience and acceptance. Experience levels will set gain targets, and noise reduction features accordingly.

The screenshot shows the 'Prescribed First Fit' settings screen. The 'Target Rule' is set to 'BAFA'. The 'Experience Level' is set to 'First Time User - Onboarding' with a sub-note: 'Focus on comfort across all input levels & noise reduction settings'. The 'Gain Level (%)' is set to '70%'. At the bottom, there is a checkbox for 'Recalculate Target Curves Only' with a sub-note: 'Check to keep existing gain curves. Best to use new target curves while performing fine tuning.'

First Time User—Onboarding

This profile is designed as a solution for patients who are new to sound amplification and need to prioritize comfort and acceptance. The prescribed gain is set to 70% of First Time User targets by default. The experience level should then be increased after an adaptation period.

The screenshot shows the 'Prescribed First Fit' settings screen. The 'Target Rule' is set to 'BAFA'. The 'Experience Level' is set to 'First Time User' with a sub-note: 'Focus on comfort for loud inputs & increased noise reduction settings'. The 'Gain Level (%)' is set to '100%'. At the bottom, there is a checkbox for 'Recalculate Target Curves Only' with a sub-note: 'Check to keep existing gain curves. Best to use new target curves while performing fine tuning.'

First Time User

This profile is designed for users who are new to hearing aids and require comfort in loud or noisy environments.

The screenshot shows the 'Prescribed First Fit' settings screen. The 'Target Rule' is set to 'BAFA'. The 'Experience Level' is set to 'Experienced - Full-time, Non-Linear' with a sub-note: 'Standard gain prescription for selected fitting rule'. The 'Gain Level (%)' is set to '100%'. At the bottom, there is a checkbox for 'Recalculate Target Curves Only' with a sub-note: 'Check to keep existing gain curves. Best to use new target curves while performing fine tuning.'

Experienced Non-Linear


This profile is for previous users or those accepting of full prescriptive targets. After an adaptation period, the patient should be set to this profile if not already done.

The screenshot shows the 'Prescribed First Fit' settings screen. The 'Target Rule' is set to 'BAFA'. The 'Experience Level' is set to 'Experienced - Full-time, Linear' with a sub-note: 'Focus on power with minimal noise reduction settings'. The 'Gain Level (%)' is set to '100%'. At the bottom, there is a checkbox for 'Recalculate Target Curves Only' with a sub-note: 'Check to keep existing gain curves. Best to use new target curves while performing fine tuning.'

Experienced Linear


This profile is designed for patients with previous experience wearing linear amplification as it provides little compression across sound input levels.

Hearing Aid Connection Flow



- 1 **Connect** the hearing aids by clicking on Connect to Beltone Solus Max in the home screen. On the discovery screen identify and assign the hearing aids from the center list by pressing the button on the hearing aid, or clicking Assign to right or left.

Note: Beltone SureFit 3 receivers will automatically select the receiver side, power and length.




- 2 Select the earpieces fitted to the hearing aids. The gray recommendation badge is provided to suggest an earpiece based on the audiogram.




- 3 To ensure a fitting without feedback, MSG calibration should be performed for every fitting. MSG calibration is an in-situ measurement of the hearing aid in the patient's ear which takes the unique parameters of the patient's hearing and ear into account. A deeper MSG curve indicates a more open fitting, thus the MSG system will be working harder.

Note: MSG calibration is required before saving for hearing aids fitted with the Microphone & Receiver-In-Ear (M&RIE) receiver.




Advanced Features and Accessories



4 Advanced features


Features which provide better access to speech are found within the Speech tab. The Comfort tab contains primarily noise reduction features as well as feedback management. The Tinnitus Breaker Pro (TBP) feature can be found in a separate tab.

Note: The hearing loss thresholds, MSG calibration, hearing aid style, and selected experience level guide the default settings for some advanced features.




5 The Device Controls menu tabs contain hearing aid hardware setup options including button configuration, beep notifications, and wireless accessory management. Within the Accessories page, you can pair with the new TV-Streamer+ and Multi-Mic+ as well as previous generation accessories. Features associated with the accessories and mobile device streaming can be adjusted here as well.

Note: Beltone wireless accessories can also be paired to the hearing aids outside of Beltone Solus Max. All paired accessories and Bluetooth® devices will appear within the Accessories tab once the hearing aids are connected.



6 For wearers without loudness sensitivity and those that have a music program, consider using Instrument MPO to avoid distortion and providing the largest possible dynamic range.

7 Click **Save** and you are done.



Additional Information

Beltone Adaptive Fitting Algorithm: The Beltone fitting rule

The BAFA (Beltone Adaptive Fitting Algorithm) target calculation is based on the basic idea of loudness normalization. To achieve loudness normalization, frequency and input level dependent amplification is applied, so that hearing aid users perceive the loudness of narrowband signals in a similar way to a person whose natural hearing dynamics are not restricted.

Hearing aid users generally tend to prefer less amplification than the basic principle of loudness normalization would predict (Smeds, 2006; Keidser & Grant, 2003). For this reason, BAFA calculates less amplification than a strict principle of loudness normalization would require.


Compared to the generic fitting rule NAL-NL2, which aims to maximize speech intelligibility, BAFA calculates 3 to 10 dB less insertion gain depending on the frequency and hearing loss.

Less listening effort with M&RIE

Each listener uses unique cues as sound travels through the environment. M&RIE gives the hearing aid wearer access to these natural cues again. The microphone positioned in the ear canal allows the hearing aid to access sound in a very individual way – reflected by the user's own ear and head anatomy. This not only makes the hearing aid sound more natural, but also provides better localization cues. Hearing aid users can better orient themselves in space and differentiate between different sound sources. The M&RIE microphone is protected by the outer ear, and as a result, wind noise is reduced by up to 22 dB compared to omnidirectional microphone transmission*.

When considering fitting M&RIE, remember the receiver is slightly larger than the LP and MP due to the additional microphone and should be fitted with a more closed fitting. The M&RIE receiver was developed with BAFA and it's therefore the recommended fitting rule for use with M&RIE.

**Schumacher (2023); Newall et al. (2023)



Display Options


Display Options contains the Vent Correction and Output View. Both displays can be valuable in a fitting.

Output View displays the Max Power Output (MPO) curve. If the MPO and gain curves are compressed together, consider increasing the MPO or changing to a higher power receiver. If adjustments were made to low frequency gain and are not audible to the patient or are not visible on real ear verification, select Vent Correction View to assist in evaluating the current earpiece selection. If the gain curves are significantly below target, consider a more occluding earpiece.



Compare across technology levels

The catalog feature has become even smarter in Beltone Solus Max 2.0. Now, you can not only look up and check all our hearing aid models, but you can also compare the differences across technology levels.

Compare hearing aid features on different technology levels for better counselling or look up the individual hearing loss for the optimal hearing aid style.



- 1 Go to “Catalog” and browse through all hearing aid models or use the filter.
- 2 Click “Right” on the selected hearing aid/technology level and “Left” on the hearing aid you would like to compare with.
- 3 Explore the different features.



US Headquarters
Manufacturer According to FDA:
Beltone
8001 E Bloomington Fwy
Bloomington, MN 55420
U.S.A.
Tel: +1 800 621 1275
Beltone.com



© 2024 GN Hearing A/S. All rights reserved. Beltone is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone and iPad and are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.