

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
3. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.
4. Social Workers: For additional information regarding standards and indicators for cultural competence, please review the NASW resource: [Standards and Indicators for Cultural Competence in Social Work Practice](#)
5. Need Help? Select the “Help” option in the member dashboard to access FAQs or contact us.

How to print Handouts

On a Mac

- Open PDF in Preview
- Click File
- Click Print
- Click dropdown menu on the right “preview”
- Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.

On a PC

- Open PDF
- Click Print
- Choose # of pages per sheet from dropdown menu
- Choose Black and White from “Color” dropdown

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

A solid red square is located in the bottom left corner of the page.

Tinnitus Demystified-A Look at Starkey's Options

Jake Hall, Au.D. Education and Training Audiologist



Hear better. Live better.

Tinnitus Demystified-A Look at Starkey's Options

Webinar FAQs

How to Ask a Question

Click Q & A, type your question,
click Send

Where to get the Handout

A PDF handout for this presentation
may be accessed through your
Pending Courses on Audiology
Online or by clicking on the link in
Chat within the classroom

Tinnitus Demystified-A Look at Starkey's Options

Webinar FAQs

How to Earn CEUs

CEUs are available with a paid Audiology Online Unlimited CEU Access membership. After completing the course, go to your Pending Courses on Audiology Online and pass the multiple-choice exam.

Assistance or Tech Support

- Call Audiology Online at 800-753-2160
- customerexperience@continued.com
- Use the Q & A window

Risks/Benefits

- The applicability and accuracy of the course content may vary by the clinical, geographical, and cultural context; you should evaluate the course accordingly.
- As healthcare is rapidly evolving, new findings may change the validity of the information presented in this course. Providers should always consult the latest research and guidelines from professional associations.
- The interventions discussed in this course may be limited in generalizability, requiring practitioners to consider specific individual and population factors.
- The course may not cover all possible risk factors, diagnostic methods, or treatment options, and complex cases may require additional considerations.
- Healthcare providers should apply cultural competence and sensitivity to ensure effective and respectful care based on their patients' diverse needs.
- Providers should always evaluate the limitations of diagnostic and screening tools, considering their potential for false results as well as any ethical implications.
- It is the responsibility of course participants to critically evaluate the evidence from multiple sources to guide professional practice.

Learner Outcomes

1

Discuss what patient controls for tinnitus are available in My Starkey app

2

Discuss how to adjust tinnitus controls in Pro Fit

3

Describe the sound therapy stimulus options available from Starkey


tin·ni·tus

["TIN-it-us", "tin-NITE-us"]

a symptom characterized by
the perception of sound in the
absence of an external stimulus



Tinnitus



One of the
most common health conditions
in the United States

American Academy of Otolaryngology
found that **one in ten Americans**
experience some form of tinnitus

The NIDCD estimates that
approximately
two million Americans
experience debilitating tinnitus

Classifying Tinnitus

Traditionally classified by whether it could be heard by the clinician



Objective

Audible to another person

Internal acoustic source

Ex. muscle spasms or vascular tumor

Detected via human ear,
stethoscope, or microphone

Represents less than 1% of cases



Subjective

Only heard by the patient


May be idiopathic or caused by
various disorders involving the
head/neck or other systems

99% of all cases

Treatment



Hear better. Live better.



Sound Therapy

Tinnitus sound therapy uses a process known as **habituation** to retrain the way the brain interprets tinnitus



Sound Therapy Options

MASKING

DISTRACTION

HABITUATION

NEUROMODULATION

Sound Therapy Devices

Sound masking devices

Hearing aids

Modified-sound/
Notched-music devices

Sound and sleep apps

Maskers in Hearing Aids

Acoustic Therapy is using sound stimulation either from amplification of external sounds or the generation of additional sounds

Tinnitus maskers are additional sounds generated that are typically set at a level below a patient's tinnitus to potentially cover it up



Masks in Hearing Aids


RESEARCH DATA

In this study, participants were provided multiple masker options to use at their discretion in order to achieve tinnitus relief

Study was comprised of 26 individuals with hearing loss

Participants completed an 8-week field trial wearing hearing aids providing acoustic therapy via three tinnitus masker options set just below minimum masking level

Tinnitus Handicap Inventory was used before starting acoustic therapy and posttreatment to measure the patients' tinnitus handicap



Source: *Journal of Speech, Language, and Hearing Research*, Vol. 64, 1415, Figure 2, April 2021.

Maskers in Hearing Aids

DATALOGGING RESULTS


Average hearing aid use was 9.4 hr (SD = 4.2 hr) during the field trial, and

Average of 34.4% of that time was spent utilizing a tinnitus masker

13.8% with white noise;

10.5% with audiogram-shaped noise;

10.0% with custom noise



Source: *Journal of Speech, Language, and Hearing Research*, Vol. 64, 1415, Figure 2, April 2021.



Multiflex Tinnitus Pro

**Next generation feature for sound
therapy tinnitus management**

Three stimulus options to
effectively mask the user's tinnitus:

- White Noise
- Audiogram-shaped
- Custom

GENESIS^{AI}

Signature
SERIES

EDGE^{AI}


AVAILABLE ON ALL STYLES

Multiflex Tinnitus Pro


Accessed from the Fly Out Menu
within Pro Fit

Use the dropdown menu to select
the stimulus type


All options offer a 10,000 Hz
bandwidth with 16 bands of
adjustability




Starkey **Pro Fit** Tinnitus Stimuli



White Noise
Pure Tone Average
All Styles and Technology Tiers



Audiogram Shaped
Hearing Loss
All Styles of the 24, 20, and 16
Technology Tiers




Custom
Custom Shaping
All Styles of the 24 and 20
Technology Tiers

Tilt Control

Low Frequency Emphasis


Click to decrease high frequency output while increasing low frequency output.



Start

High Frequency Emphasis

Click to increase high frequency output while decreasing low frequency output.




Start

Tilt Control

Low Frequency Emphasis


Click to decrease high frequency output while increasing low frequency output.



5 clicks

High Frequency Emphasis

Click to increase high frequency output while decreasing low frequency output.




5 clicks

———— = Original response


Starkey **Pro Fit**

Customizing Noise Colors



Tilt

< >




Tinnitus Fine Tuning

Tinnitus Fine Tuning						
Low	Mid	High	Expand All	Low	Mid	High
150 - 1.3k	1.9k - 3.1k	3.8k - 9.4k	ALL	150 - 1.3k	1.9k - 3.1k	3.8k - 9.4k
48	16	-6		48	16	-6

Modulation

Modulation is available for all the Multiflex Tinnitus stimulus options



Starkey Pro Fit


White Noise Tinnitus

Signal level based on the PTA and comprised of equal energy distribution across the frequencies

Flat response may result in some segments of the stimulus being too loud or too soft

Vent interactions are accounted for

Stimulus is adjustable




Starkey **Pro Fit**

Audiogram Shaped Tinnitus

Automatic shaping of the stimulus based on the hearing loss provides a more tailored masking stimulus

Vent interactions are accounted for

Stimulus is adjustable




Starkey Pro Fit

Custom Tinnitus

Custom shaping of the stimulus provides a personalized stimulus for a greater degree of masking precision

Vent interactions are accounted for

Stimulus is adjustable



Custom

Personalized shaping of the Multiflex Tinnitus stimulus

Multiflex Tinnitus Pro Tool is a research-based approach to improving the efficiency of the masking signal

Minimum Detection Levels (MDLs) and Minimum Masking Levels (MMLs) provide personalized shaping of the tinnitus stimulus



Efficient and Effective Masking

Stimulus Personalization Testing in Pro Fit typically only takes a few minutes and is less time-consuming than tinnitus measurements like pitch-matching

MDL and MML measurements can be used to rapidly design a stimulus that is both efficient and effective at masking tinnitus

To efficiently mask tinnitus, only those frequency bands that are effective at masking tinnitus should be included in the stimulus, and the level of the stimulus in these bands should be no greater than needed to achieve masking


Fournier et al., 2018

Starkey Pro Fit

Custom

Tinnitus

Stimulus Personalization



Minimum Detection Level Test

- To enable presentation of the masker, click the Stimulus button
- Presentation of the stimulus is constant when enabled
- First, determine stimulus Minimum Detection Level (MDL) for each of three frequency groups
- Click the arrow to move to the next frequency group
- Once MDL is measured for all three frequency groups, select Next to advance to the Minimum Masking Level (MML) test

[? more details](#)



Next

Cancel

Stimulus Personalization



Minimum Masking Level Test

- To enable presentation of the masker, click the Stimulus button
- Presentation of the stimulus is constant when enabled
- First, determine stimulus Minimum Detection Level (MDL) for each of three frequency groups
- Click the arrow to move to the next frequency group
- Once MDL is measured for all three frequency groups, select Next to advance to the Minimum Masking Level (MML) test

[? more details](#)



Could Not Test

Could Not Test

Apply Binaural

Cancel

Multiflex Tinnitus Pro Tool Display

Active frequency region

Instructions



Stimulus Controls

Stimulus Personalization

Right

Binaural

Left



Minimum Detection Level Test

- To enable presentation of the masker, click the Stimulus button
- Presentation of the stimulus is constant when enabled
- First, determine stimulus Minimum Detection Level (MDL) for each of three frequency groups
- Click the arrow to move to the next frequency group
- Once MDL is measured for all three frequency groups, select Next to advance to the Minimum Masking Level (MML) test

[? more details](#)




Next

Cancel


Frequency Groups

Frequency groupings are consistent with the low, mid and high band groupings in Inspire




Low

150 Hz to 1.3 kHz



Mid

1.9 kHz to 3.1 kHz





High

3.8 kHz to 9.4 kHz

Minimum Detection Level (MDL) Test

Stimulus Personalization


Right Binaural Left



Minimum Detection Level Test

- To enable presentation of the masker, click the Stimulus button
- Presentation of the stimulus is constant when enabled
- First, determine stimulus Minimum Detection Level (MDL) for each of three frequency groups
- Click the arrow to move to the next frequency group
- Once MDL is measured for all three frequency groups, select Next to advance to the Minimum Masking Level (MML) test

[? more details](#)



GOAL

To identify the softest detectable level of the Multiflex Tinnitus stimulus

- 1 Click on the 'Stimulus' button to turn it ON
- 2 Hearing aid microphones will mute automatically when the stimulus is ON
- 3 Use the up and down arrows to increase or decrease the level of the stimulus in the frequency group displayed to **identify when the stimulus is just detectable**
- 4 Use the right and left arrows to make measurements in the other frequency groups


Minimum Masking Level (MML) Test

GOAL

To identify the softest level of the Multiflex Tinnitus stimulus required to just notice or mask the tinnitus depending on the management approach of choice

Display and controls are the same for this test

Increase the stimulus level to determine ***the softest stimulus level required to just detect or mask the tinnitus*** for each of the three frequency groups



Tinnitus Stimulus

User Controls

● Short Press

Tinnitus Level ✕

☒ Ear to Ear →

● Push & Hold

Program ✕

☒ Ear to Ear →

● Double Tap

Edge Mode+ ✕

☒ Ear to Ear →

● Short Press

Volume ✕

← ☒ Ear to Ear

● Push & Hold

Dissonance

☐ Ear to Ear

● Double Tap

Smart Assistant ✕

☐ Ear to Ear

Output Range and Step Size

User Control Settings

Manage Phone Calls

VC Range and Step Size

Tinnitus Output Range and Step Size

Tinnitus Output Range: -6dB to 6dB

Step Size

2 dB

4 dB

-6dB

6dB

Tinnitus Output Range: -6dB to 6dB

Step Size

2 dB

4 dB

-6dB

6dB

Additional Configuration Options


Pro Fit Configuration

Multicurve View:

Multiple “Personal” Programs

Varying amounts of gain for fluctuating hearing loss

Tinnitus Masker- Audiogram shaped with slow modulation for all four programs



Multiflex Tinnitus Pro


Personalized Sound Therapy for Tinnitus Management



WHITE NOISE

Signal level based on PTA and comprised of equal energy distribution across frequencies


All Styles and Technology Tiers



AUDIOGRAM-SHAPED NOISE

Automatic shaping of the tinnitus stimulus based on the hearing loss

All Styles of the 24, 20 and 16 Technology Tiers




CUSTOM NOISE

Custom shaping of the tinnitus stimulus based on MDLs and MMLs

All Styles of the 24 & 20 Technology Tiers


Multiflex Tinnitus Pro

Low-frequency weighted options




OCEAN NOISE

Mimics the sound of the ocean for an alternative masker to traditional options



PINK NOISE

Comprised of less high-frequency energy compared to White Noise




RED NOISE

Comprised of less high-frequency energy compared to Pink Noise

All Styles and Technology Tiers of Edge AI and Signature Series


Multiflex Tinnitus Pro

High-frequency weighted options



BLUE NOISE

Comprised of less low-frequency energy compared to White Noise



PURPLE NOISE

Comprised of less low-frequency energy compared to Blue Noise

All Styles and Technology Tiers of Edge AI and Signature Series




Tech Tier and Style Differences


Edge AI & Signature Series

Genesis AI

White Noise	All styles and tiers	All styles and tiers
Audiogram-shaped	All styles and tiers	All styles; 24, 20, & 16
Custom Masking	All styles; 24 & 20	All styles; 24 & 20
Color Masking	All styles and tiers	N/A
Ocean	All styles and tiers	N/A

My Starkey Controls





Relief sounds and
education for patients
who have tinnitus



Starkey Relax

Set Sound for Masking

Set Modulation, Volume, and Timer

The quickest way to ruin a customer experience in healthcare is to treat everyone the same

Patients don't want to feel like just another number

They crave personalized service that helps them find the right solutions

Hyken, S.

Personalized Customer Experience Increases Revenue And Loyalty.

Tinnitus Resources

American Tinnitus Association

www.ata.org

Tinnitus Research Initiative

www.tinnitusresearch.net

**National Center for Rehabilitative Auditory
Research (NCRAR)**

[www.ncrar.research.va.gov/Education/Documents
/TinnitusDocuments](http://www.ncrar.research.va.gov/Education/Documents/TinnitusDocuments)

American Academy of Audiology

www.audiology.org

Tinnitus Retraining Therapy

www.tinnitus-pjj.com

**University of Iowa Tinnitus & Hyperacusis
Research**

[www.medicine.uiowa.edu/oto/research/tinnitus-
and-hyperacusis](http://www.medicine.uiowa.edu/oto/research/tinnitus-and-hyperacusis)

Resources

StarkeyPro.com

- Multiflex Tinnitus Pro QuickTIP
- Multiflex Tinnitus Technology Handbook
- Tinnitus Handicap Inventory Questionnaire
- Tinnitus Consumer In Office Poster
- Patient Education Brochure: Tinnitus
- Tinnitus Technology Options Brochure

TINNITUS HANDICAP INVENTORY

Patient Name: _____ Date: _____

INSTRUCTIONS: The purpose of this questionnaire is to identify difficulties that you may be experiencing because of your tinnitus. Please answer every question. Please do not skip any questions.

1. Because of your tinnitus, is it difficult for you to concentrate?	Yes	Sometimes	No
2. Does the loudness of your tinnitus make it difficult for you to hear people?	Yes	Sometimes	No
3. Does your tinnitus make you angry?	Yes	Sometimes	No
4. Does your tinnitus make you feel confused?	Yes	Sometimes	No
5. Because of your tinnitus, do you feel desperate?	Yes	Sometimes	No
6. Do you complain a great deal about your tinnitus?	Yes	Sometimes	No
7. Because of your tinnitus, do you have trouble falling to sleep at night?	Yes	Sometimes	No
8. Do you feel as though you cannot escape your tinnitus?	Yes	Sometimes	No
9. Does your tinnitus interfere with your ability to enjoy your social activities, such as going out to dinner, to the movies?	Yes	Sometimes	No
10. Because of your tinnitus, do you feel frustrated?	Yes	Sometimes	No
11. Because of your tinnitus, do you feel that you have a terrible disease?	Yes	Sometimes	No
12. Does your tinnitus make it difficult for you to enjoy life?	Yes	Sometimes	No
13. Does your tinnitus interfere with your job or household responsibilities?	Yes	Sometimes	No
14. Because of your tinnitus, do you find that you are often irritable?	Yes	Sometimes	No
15. Because of your tinnitus, is it difficult for you to read?	Yes	Sometimes	No
16. Does your tinnitus make you upset?	Yes	Sometimes	No
17. Do you feel that your tinnitus problem has placed stress on your relationships with members of your family and friends?	Yes	Sometimes	No
18. Do you find it difficult to focus your attention away from your tinnitus and on other things?	Yes	Sometimes	No
19. Do you feel that you have no control over your tinnitus?	Yes	Sometimes	No
20. Because of your tinnitus, do you often feel tired?	Yes	Sometimes	No
21. Because of your tinnitus, do you feel depressed?	Yes	Sometimes	No
22. Does your tinnitus make you feel anxious?	Yes	Sometimes	No
23. Do you feel that you can no longer cope with your tinnitus?	Yes	Sometimes	No
24. Does your tinnitus get worse when you are under stress?	Yes	Sometimes	No
25. Does your tinnitus make you feel insecure?	Yes	Sometimes	No

FOR CLINICIAN USE ONLY

Total Per Column	x4	x2	x0
Total Score			

To receive the score please refer to the Tinnitus Handicap Inventory Scale shown on the inside back.



Starkey Resources

StarkeyPro.com

JSLHR

Research Note

Changes in Heart Rate Variability Following Acoustic Therapy in Individuals With Tinnitus

Paul Reinhart,^a Kendra Griffin,^a and Christophe Micheyl^b

Purpose: The aim of this study was to investigate changes in autonomic function, as measured by heart rate variability, in individuals with tinnitus following acoustic therapy. Implemented using tinnitus maskers presented via hearing aids.

Method: Twenty-six individuals with tinnitus and hearing impairment completed an 8-week field trial wearing hearing aids providing acoustic therapy via three tinnitus masker options set just below minimum masking level. Tinnitus handicap was measured using the Tinnitus Handicap Inventory at baseline (before starting acoustic therapy) and posttreatment (at end of 8-week trial). Resting heart rate and

heart rate variability were measured using electrocardiography at baseline and posttreatment.

Results: There was a significant decrease in tinnitus handicap posttreatment compared to baseline. There was no change in heart rate, but there was a significant increase in heart rate variability posttreatment compared to baseline.

Conclusions: Acoustic therapy using tinnitus maskers delivered via hearing aids provided tinnitus relief and produced a concurrent increase in heart rate variability, suggesting a decrease in stress. Heart rate variability is a potential biomarker for tracking efficacy of acoustic therapy; however, further research is required.

Tinnitus refers to auditory perceptions that do not have a corresponding acoustic source in the external environment. Estimates of the prevalence of tinnitus vary widely, ranging from 3% to 43%, depending on the population sampled, and what type of tinnitus is considered (e.g., fleeting or chronic, bothersome or not; McCannack et al., 2016). Chronic tinnitus can interfere with the ability to fall asleep or concentrate, and often induces stress, anxiety, and depressive feelings, in turn reducing quality of life (Wack et al., 2016). It has been hypothesized that tinnitus-related distress is related to an imbalance of sympathetic and parasympathetic branches of the autonomic nervous system, which modulate active and rest states, respectively (Jastreboff & Jastreboff, 2000). Consistent with this hypothesis, individuals with tinnitus exhibit decreased autonomic function (increased sympathetic activity and suppressed parasympathetic activity) indicative of being in a constant elevated stress state (Bartov et al., 1999;

van Der Loos et al., 2011; Vanneste & De Ridder, 2013; Ylikoski et al., 2017).

Heart rate variability (HRV) is a common, noninvasive measure reflecting sympathetic and parasympathetic branches of the autonomic nervous system (Bett et al., 2017; O'Brien et al., 1986; Task Force of the European Society, 1996). The heart, which is regulated by the autonomic nervous system (Glick et al., 1965; Valsarta et al., 2018), contracts with slight beat-to-beat irregularity; the quantification of this inherent irregularity is referred to as HRV (Acharya et al., 2006). HRV is commonly used as an autonomic biomarker of stress and health, with lower HRV associated with increased stress and negative health outcomes (Kemp & Quintana, 2013; Thayer et al., 2012).

Acoustic therapy, also known as "sound stimulation" or "sound therapy," constitutes an important aspect of comprehensive tinnitus management strategies used in audiology today (Henry et al., 2018; Jastreboff & Jastreboff, 2000; Schud et al., 2018; Tunkel et al., 2014). Acoustic therapy involves sound stimulation resulting, either or both, from the amplification of external sounds, or the generation of additional sound, commonly referred to as a "tinnitus masker" or "tinnitus stimulus." Tinnitus maskers are commonly set to a level below their tinnitus, so that their tinnitus is partially covered up (Jastreboff & Jastreboff, 2000). Because hearing loss is a common comorbidity for tinnitus

new to acoustic therapy via hearing aid tinnitus maskers. Individuals with hyperacusis, as indicated by a comorbidity score > 4 on the Social Tolerance section of the Tinnitus and Hearing Survey, were excluded from the study (Henry et al., 2018). Participants were in good general health at time of testing and did not report any history of cardiovascular impairment. Participants completed an informed consent process and were compensated for their time.

Acoustic Therapy

Acoustic therapy was implemented using tinnitus maskers (aka "tinnitus stimulus") presented via hearing aids. Participants were fit with Starkey 24-channel hearing aids—five participants received behind-the-ear devices with custom earmolds, and 21 participants received receiver-in-canal devices with either open (3/21), closed (16/21), or power domes (2/21), depending on audiometric configuration. All fittings were bilateral, except for one participant with unilateral hearing loss. Hearing aid gain was calculated using the NAL-NL2 prescriptive formula (Kiddner et al., 2011) and was verified by matching the real-ear aided response to a ± 2 -dB target using the International Speech Test Signal (Hohle et al., 2010) presented at soft (50 dB SPL), average (65 dB SPL), and loud (80 dB SPL) speech levels. In addition, user controls were enabled so that participants could adjust the hearing aid volume in 4-dB increments via remote control accessory.

There is no one-size-fits-all approach to tinnitus acoustic therapy, with different maskers providing varying relief across individuals (e.g., McNeill et al., 2012; Scherette et al., 2010). As such, the study hearing aids were programmed with three different tinnitus maskers available as part of Starkey's Multiflex Tinnitus Pro technology:

1. White noise—a noise with nominally flat spectrum (i.e., constant sound level per band across frequencies).
2. Audiogram-shaped noise—a noise spectrally shaped based on individual participants' audiogram to provide approximately equal audibility across all frequencies.
3. Custom noise—a noise based on psychophysical measurements of minimum masking level (MML), which were used to generate a customized masker spectrum primarily containing energy at frequencies heard where masking is most effective (Fournier et al., 2018; Vernon & Merkle, 2003).

The reason for providing multiple masker options was so that participants could utilize the acoustic therapy option that provides them the most tinnitus relief (Konttinen et al., 2020). Any potential changes in HRV following acoustic therapy were agnostic to the specific tinnitus masker commonly set to a level below their tinnitus, so that their tinnitus relief was achieved. Participants could listen to the different maskers by changing memories via hearing aid rocker switch or remote control accessory. The level of each masker was initially set to just below MML (i.e., the level at which their tinnitus was completely masked), consistent with tinnitus

Results


Demographic results extracted at the end of the study indicated that average hearing aid use was 9.4 h ($SD = 2.3$ hr) during the field trial, and an average of 34.4% of test time was spent utilizing a tinnitus masker (13.8% with this noise; 10.2% with audiogram-shaped noise; 10.0% with custom noise).

THI, heart rate, and HRV data were analyzed using paired-samples t tests to compare data at baseline (beginning of study) versus posttreatment (end of 8-week trial). There was a significant decrease in THI posttreatment ($M = 8.3$, $SD = 14.5$) compared to baseline ($M = 48.8$, $SD = 8.7$, $t(24) = 6.046$, $p < .001$, Cohen's $d = 1.385$ (see Figure 2). There was no significant change in heart rate between baseline ($M = 67.4$, $SD = 11.7$) and posttreatment ($M = 67.4$, $SD = 9.3$, $t(24) = 1.022$, $p = .316$ (see Figure 3, left panel). There was a significant increase in HRV posttreatment ($M = 50.2$, $SD = 57.3$) compared to baseline ($M = 33.5$, $SD = 35.5$, $t(24) = -2.142$, $p < .05$, Cohen's $d = 0.430$ (see Figure 3, right panel). A Pearson correlation was conducted to examine the relationship between change in HRV baseline HRV—posttreatment HRV and change in THI baseline THI—posttreatment THI. This test failed to reach significance, $r(23) = .124$, $p = .322$.

Discussion


In this study, participants exhibited a decrease in tinnitus handicap and an increase in HRV following an 8-week acoustic therapy field trial. This finding was consistent with our prediction and supported the hypothesis that acoustic therapy induces an increase in parasympathetic activity in order to restore autonomic sympathetic-parasympathetic balance in individuals with chronic, bothersome tinnitus. Previous tinnitus research has found a correlation between HRV and activity in the dorsal and subgenual anterior cingulate cortex, suggesting that those brain regions could be

Figure 2. Results of Tinnitus Handicap Inventory showing a significant reduction ($p < .001$) in tinnitus handicap posttreatment ($M = 8.3$) compared to baseline ($M = 48.8$) (beginning of study); error bars represent ± 1 standard error.



Reinhart et al. / Heart Rate Variability and Tinnitus 1415

Figure 3. (Beginning of study) and posttreatment (end of 8-week trial). Left: Resting heart rate variability (in non-normalized units of successive differences).




control by being able to choose among three masker options. This would suggest that tinnitus patients might benefit most when provided multiple masker options; however, this requires further research.

HRV has had limited use in audiology and has primarily been used as a measure of acute listening-related stress occurring during challenging listening speech perception, such as speech in noise (Dumas et al., 2012; Mackenzie et al., 2015; Scemman & Sim, 2015) or time-compressed speech (Mackenzie & Calderon-Morales, 2016). While validated self-report measures, such as THI (Newman et al., 1998), Tinnitus Reaction Questionnaire (Wilson et al., 1991), or Tinnitus Functional Index (Mackenzie et al., 2012), and THI baseline THI—posttreatment THI, this test failed to reach significance, $r(23) = .124$, $p = .322$.

Given the potential relationship between tinnitus handicap and HRV, patients with tinnitus may be able to benefit from HRV biofeedback training (Leherer & Gevirtz, 2014). This training teaches individuals to maximize HRV resulting from respiratory artery arrhythmia through controlled breathing (Leherer et al., 2011). HRV biofeedback training has been shown to provide clinical benefits in a range of conditions, such as pain (Hullman et al., 2011) and anxiety (Reiner, 2008). It is possible that individuals with tinnitus could also benefit from HRV biofeedback training, as other modulation and relaxation therapies have shown benefit in individuals with chronic, bothersome tinnitus (An et al., 2017).

In conclusion, this study found a significant increase in HRV, corresponding to a significant decrease in tinnitus



Starkey
Hearing Technologies

^aStarkey Hearing Technologies, Eden Prairie, MN
^bStarkey Tinnitus Center
 Correspondence to: Paul Reinhart (paul_reinhart@starkey.com)
 Editor-in-Chief: Frederick Frenk (frenk@starkey.com)
 Editor: Alexander L. Frenk
 Received October 9, 2020
 Revision received January 9, 2021
 Accepted January 11, 2021
 https://doi.org/10.1044/2021.JSLHR-20-00796

Thank you



Hear better. Live better.