


Voice Aligned Compression (VAC+)

Oticon's proprietary fitting rationale



Virginia Ramachandran, AuD, PhD
Vice President, Audiology, Oticon, Inc.

oticon
life-changing technology

Virginia Ramachandran, AuD, PhD
Vice-President, Audiology, Oticon, Inc.
vira@Oticon.com



Speaker Disclosure

Relevant financial relationships:

- Employee of Oticon, Inc. and receives a salary
- Author & Associate Editor, Plural Publishing, Inc.

Non-financial relationships:

- Past-president American Academy of Audiology

Learner Outcomes

1. After this course, participants will be able to describe two key benefits of the Voice Aligned Compression (VAC+) proprietary fitting rationale.
2. After this course, participants will be able to list the options available to verify Oticon's proprietary fitting rationale when completing real ear measures, including use of the SII.
3. After this course, participants will be able to describe the results of a pilot study examining the audibility provided by VAC+ and measures of subjective experience.

Risks / benefits

- The applicability and accuracy of the course content may vary by the clinical, geographical, and cultural context; you should evaluate the course accordingly.
- As healthcare is rapidly evolving, new findings may change the validity of the information presented in this course. Providers should always consult the latest research and guidelines from professional associations.
- The interventions discussed in this course may be limited in generalizability, requiring practitioners to consider specific individual and population factors.
- The course may not cover all possible risk factors, diagnostic methods, or treatment options, and complex cases may require additional considerations.
- Healthcare providers should apply cultural competence and sensitivity to ensure effective and respectful care based on their patients' diverse needs.
- Providers should always evaluate the limitations of diagnostic and screening tools, considering their potential for false results as well as any ethical implications.
- It is the responsibility of course participants to critically evaluate the evidence from multiple sources to guide professional practice.

Agenda

1. Introduction
2. Assumptions underlying rationales
3. The evolution of VAC+
4. VAC+ in Genie2
5. Verification of VAC+
6. Pilot study results of VAC+ research
7. Conclusions

What is the **main goal**
of a **hearing aid fitting**?

**To fit hearing aids that meet your client's
unique needs and preferences**



oticon



Steps for successful communication



Adapted from Kiessling et al., 2003

Fitting rationales

Algorithms that guide the prescription of gain




oticon

Fitting rationales and compression

Provide gain targets based upon the client's characteristics

Example of
VAC+ targets



“In general, prescriptive methods are based on the assumption that using the patient’s pure-tone thresholds, and perhaps some supra-threshold measures, hearing aid gain and/or output can be *prescribed* that will result in the best fitting **for the average patient** for the desired attributes (e.g. intelligibility, sound quality, listening comfort)...With modern hearing aids the intent is also to accomplish these goals for a wide range of sound input levels and listening environments. There are a number of variables to consider, and no single procedure is going to get it perfect for all patients all the time”.

-Ricketts, T. A., Bentler, R., & Mueller, H. G. (2019). Essentials of Modern Hearing Aids: Selection, Fitting, & Verification. Plural Publishing, Inc.


A possible misconception



So, if I perfectly
**match prescription
targets**, my patient
will be happy with the
sound of their
hearing aids?

“In general, prescriptive methods are based on the assumption that using the patient’s pure-tone thresholds, and perhaps some supra-threshold measures, hearing aid gain and/or output can be *prescribed* that will result in the best fitting for the average patient for the desired attributes (e.g. intelligibility, sound quality, listening comfort). With modern hearing aids the intent is also to accomplish these goals for a wide range of sound input levels and listening environments. There are a number of variables to consider, and no single procedure is going to get it perfect for all patients all the time”.

-Ricketts, T. A., Bentler, R., & Mueller, H. G. (2019). Essentials of Modern Hearing Aids: Selection, Fitting, & Verification. Plural Publishing, Inc.





precision care

“Most medical treatments are designed for the “average patient” as a one-size-fits-all-approach, which may be successful for some patients but not for others...The goal of precision medicine is to target the right treatments to the right patients at the right time.”

<https://www.fda.gov/medical-devices/in-vitro-diagnostics/precision-medicine>



<https://www.library.ucdavis.edu/guide/ebp-resources/>

CONSUMER reports

Consumers Union test results

September 1950

WINDOW SHADES
TISSUES FACIAL & TOILET
EVAPORATED MILK
THE NEW KAISER

for complete contents see third page

35¢

DEAFNESS

AND HEARING AIDS



Again, an “audiogram”—that is, the curve showing relative loss of hearing at various frequencies—is drawn by determining the *least loud sound* of a particular pitch that you can hear. This is the “threshold of audibility.” But research has shown that a threshold audiogram may be a wholly misleading guide to the way you hear sounds louder than threshold sounds. Since the earphone delivers sounds well above threshold levels, fitting a hearing aid to your threshold audiogram was a wholly illogical procedure.

The whole elaborate theory of “selective amplification” was finally disposed of in 1947, when Dr. Hallowell Davis and his associates published their findings. They had started their Federal research project at Harvard with the usual “personalized fitting” theory, and had constructed a “master hearing aid” which could be set for a wide range of frequency response curves. They then tried various curves on a group of trained hard-of-hearing listeners, and gauged the results by means of carefully designed tests of speech intelligibility. Their finding was that two types of curves were suitable to practically all of the hard-of-hearing subjects tested.

The “selective amplification” theory is still being pushed by some companies, but CU believes that it is on its last legs. It has outlived its usefulness even to those companies which have profited most from it over the past two

HEARING AIDS

AN EXPERIMENTAL STUDY

of
DESIGN OBJECTIVES

BY

HALLOWELL DAVIS
Project Supervisor

S. S. STEVENS
Director, Psycho-Acoustic Laboratory

R. H. NICHOLS, JR.
Associate Director, Electro-Acoustic Laboratory

C. V. HUDGINS G. E. PETERSON
R. J. MARQUIS D. A. ROSS



The validation of the laboratory studies by large-scale clinical testing at one of the Army Aural Rehabilitation Hospitals, which had originally been planned, has not been undertaken. The termination of the wartime activities of the Office of Scientific Research and Development is partly responsible for this curtailment, but the unexpectedly simple and definite result of the laboratory experiments makes the large-scale validation seem less vital than was anticipated. The results also seem to show that individual testing with the Master Hearing Aid is unnecessary for the satisfactory selection of hearing aids for the majority of patients (1),¹ so that in an unexpected way the second objective has also been largely attained. We believe, in other words, that we have disproved the fundamental assumption of the desirability of "selective amplification" based on the characteristics of the individual's audiogram.

	Subject	Age	HL Type
1	RA	44	Mixed
2	FB	49	Mixed (presumed otosclerosis)
3	HB	53	Mixed (presumed otosclerosis)
4	BD	30	Mixed
5	BL	70	Mixed
6	DL	43	Mixed (presumed otosclerosis)
7	HM	40	Mixed (presumed otosclerosis)
8	RR	46	Mixed
9	JH	17	SNHL
10	IS	23	SNHL
11	RW	48	SNHL
12	PP	25	SNHL (? - only 1 BC threshold shown)
13	WW	23	SNHL (probably - no BC responses)
14	MD	?	Unknown (no BC performed)
15	TR	?	Unknown (no BC performed)
16	MC	52	Unknown (no audiogram - presumed SNHL)
17	FS	64	Unknown (no audiogram - presumed SNHL)
18	RC	37	Unknown (no audiogram - presumed conductive)


- Normal hearing function relies on 30+ different structures/functions that must work in a precisely coordinated sequence

Outer Ear	Middle Ear	Inner Ear (Cochlea)	Auditory Nerve and Spiral Ganglion	Other Key Interfaces
Pinna (Auricle) External auditory canal (meatus)	Tympanic membrane (eardrum) Malleus Incus Stapes Eustachian tube Middle ear muscles (tensor tympani and stapedius) Middle ear ligaments and joints (e.g., incudomalleolar and incudostapedial joints)	Oval window Round window Scala vestibuli Scala tympani Scala media (cochlear duct) Basilar membrane Tectorial membrane Organ of Corti Inner hair cells Outer hair cells Supporting cells (e.g., Deiters', Hensen's, Claudius' cells) Stria vascularis Spiral ligament Endolymph and perilymph fluids	Spiral ganglion neurons Peripheral processes of spiral ganglion neurons (connect to hair cells) Central axons of spiral ganglion neurons (form auditory nerve) Cochlear nerve (branch of CN VIII) Myelin sheaths (Schwann cells)	Synapses between hair cells and auditory nerve fibers Blood supply to cochlea and auditory nerve (e.g., labyrinthine artery)


Known, well-documented causes of hearing loss



Normal Hearing Loss



Dobie, 2015






A photograph of a man and a woman in a professional setting. The man, on the left, has grey hair and is wearing a blue and white checkered button-down shirt. He is looking towards the woman. The woman, on the right, has dark hair tied in a bun, wears glasses, and a grey blazer over a green top. They are both looking at each other in conversation. The background features vertical blinds on the left and a desk lamp on the right. A semi-transparent blue banner is at the bottom.

The importance of history & context

Case Example

- 30-year-old female
- Reports normal hearing until after experiencing hearing loss and tinnitus after a concert
- How would assumptions and average fitting data affect hearing aid fitting?





*“An approximate answer to the right
problem is worth a good deal more
than an exact answer to an
approximate problem”*

– John Tukey

precision care

"Most medical treatments are designed for the "average patient" as a one-size-fits-all-approach, which may be successful for some patients but not for others...The goal of precision medicine is to target the right treatments to the right patients at the right time."

<https://www.fda.gov/medical-devices/in-vitro-diagnostics/precision-medicine>

Eriksholm Research Center

<https://www.eriksholm.com/publication>






Background to fitting rationales

The challenge of prescribing gain

A wide range of input levels and listening environments

Restore
audibility





The challenge of prescribing gain

A wide range of input levels and listening environments

Restore
audibility

Maximize **speech**
intelligibility




Avoid
overamplifying

Protect
residual **hearing**


Dynamic range

Normal hearing



Dynamic range

Sensorineural hearing loss



The challenge of prescribing gain

A wide range of input levels and listening environments

Restore
audibility

Maximize **speech**
intelligibility

Preserve
sound quality



Avoid
overamplifying

Protect
residual **hearing**

Ensure
comfort

Which algorithm should I use?

It depends which
assumptions you are
willing to make.



What are some assumptions we might make?

- All speech is valuable
- Restore audibility as much as possible
- Prescriptive approach
- Correct for threshold loss
- Measurable hearing is useable hearing
- The more bandwidth the better
- Make the full range of inputs fit
- Both ears contribute equally
- Closed ear canal
- Listening in quiet
- Speech is the signal of interest



Which rationale is best for my client?

Each has a different philosophy

Generic

NAL-NL2 /
DSL v 5.0

Proprietary

VAC+

oticon



Which rationale is best for my client?

Each has a different philosophy

Generic

NAL-NL2 /
DSL v 5.0

Proprietary

VAC+

Standardized
approach

Length of time
between revisions

Built and updated
to support innovations

Historically, **limited**
information released

oticon

Generic rationale

NAL-NL2



Philosophy

Maximize speech
intelligibility



Latest update

2011



Personalization

Experience



Targeted rationales

Binaural

Ensure overall
loudness comfortable

Gender

Open fitting

Age (child / adult)

Language

Generic rationale

DSL v 5.0



Philosophy

Audibility of speech
at all frequencies

Normalize loudness
at each frequency



Latest update

2005



Personalization

Infant / child / adult

Listening
environments



Targeted rationales

Binaural

Conductive
hearing losses

Oticon's proprietary rationale

VAC+



Philosophy

Natural sound

Audibility

Comfort



Latest update

2024



Personalization

Soft sounds

Brightness



Targeted rationales

Open fitting


Complex
audiograms

The background image is a dark, futuristic tunnel. It features a series of glowing blue and orange lines that curve and recede into the distance, creating a sense of depth and motion. The lines appear to be part of a high-speed train or a digital data path. The overall atmosphere is mysterious and technological.

Evolution

The evolution of VAC+

Commitment to optimization




Ghamkhar, 2025.
Oticon Whitepaper

oticon

The evolution of VAC+

Commitment to optimization



Ghamkhar, 2025.
Oticon Whitepaper

oticon

Validation of VAC+

Continuous and rigorous internal testing

Small scale development tests



To assess **performance** with
hardware and **signal processing**

Larger scale development tests



To validate **audiological**
performance



VAC+ in practice


VAC+ in Oticon Genie 2

The **default** fitting rationale for adults

Available for **all premium** hearing aids

oticon

Personalization



51 Voice Aligned Compression (VAC+)

Personalization Parameters

Client data

Experience (right ear)	Experience (left ear)	Client language	Gender	Age
Long-term ▼	Long-term ▼	English ▼	Female ▼	60–69 years ▼

Personalization via ACT

Client

Welcome

Prepare

Fit

Finish

Client data

Experience (right ear)
None

Experience (left ear)
None

Client language
None

Gender
Unspecified

Age
60–69 years

Personalization method

Select only one.

Audible Contrast Threshold (ACT)
Non-language based test of speech-in-noise ability. The fitting software personalizes help-in-noise based on results.

Personalization questionnaire
Personalizes help-in-noise based on the client's subjective preferences.

The client has **mild** problems with understanding **speech in noise**

Mild (6 dB nCL)

Normal

Mild

Moderate

Severe

-4

0

4

7

10

16 dB nCL

How to measure ACT?

ACT is measured using an audiometer and takes 2-3 minutes. [Check here for more details about ACT.](#)

Last change: 12/23/2025

Edit value

Personalization Questionnaire

Personalization method

Select only one.

Audible Contrast Threshold (ACT)

Non-language based test of speech-in-noise ability. The fitting software personalizes help-in-noise based on results.

Personalization questionnaire

Personalizes help-in-noise based on the client's subjective preferences.

 Calculations will now be based on the questionnaire, not ACT



Adjust the volume

When using sound demos, begin by playing the reference sound demo and adjust the loudspeakers so that it is comfortable for your client.

When I wear my hearing instruments:

I am easily disturbed by sounds around me

☐ Yes

☐ Only sometimes

I prefer to hear the full range of sounds in the listening environment

 ☐ Yes

 ☐ No

I prefer sound to be

 ☐ Sharp and distinct

 ☐ Soft and round

My everyday listening environments are

 ☐ Variable

 ☐ Less variable

I prefer a more comfortable sound even if it takes away the softer details in the sound

 ☐ Yes

 ☐ No

Recalculate settings

New setting recommendations will be calculated because:

- ACT was measured

 Fine-tuning will not be affected by the changes



Recalculate new help-in-noise settings



 **Recalculate**

Keep existing settings

Keep settings


Cancel






Graph Types


P1: General (VAC+)




P1: General (VAC+)



Graph Settings



Changing prescriptive algorithms



Additional & Special Purpose Programs

×

Add program

Main programs

 General

 MyMusic

Special purpose programs

 Speech in Noise

 Comfort

 Lecture

 Phone Mic

Telecoil program

 T-coil

 T-coil & Mic

 Phone T-coil

 Phone Mic & T

MyMusic

This program is optimized to provide the best possible music experience.

Based on
Target

Internal note

Client information ⓘ


Program name

Additional information

Save program


Cancel

Gain Preferences



Personalized for your client

Brightness and soft sound perception



oticon

Client

Welcome

Prepare

Fit

Finish

Client

age: 64

Fit

P1 General VAC+

Fine-tuning

Feedback analyser

REM tools

MoreSound Intelligence

Program manager

Acoustics

Adaptation manager

Automatics

Fitting assistant

Tinnitus

Speech Rescue

Data logging

In-situ Audiometry

Intent 1 miniRITE

Bass dome, double

85

Disconnect

Save and exit

Intent 1 miniRITE

Bass dome, double

85

P1: General (VAC+)

Insertion gain

Graph settings

P1: General (VAC+)

ANSI S3.5 (ICRA Stationary)

Target 65

Simulated 65

Reserve gain 10 dB

dB

125 250 500 1k 2k 4k 8k 10k Hz

Gain controls

Sound controls

SuddenSound Stabilizer

Brightness perception

Fuller

Brighter

Soft sound perception

Comfort

Detail

Brightness perception

Fuller

Brighter


Soft sound perception

Comfort

Detail

Verifying VAC+ targets

Real-ear measurements



oticon

Important Points About REM


- Performing real ear measurement is not the same thing as matching targets.
- Real-ear is about measuring what is happening in the ear canal.
- It CAN be used to match targets
- Regardless of which algorithm you choose to use and whether or not your clinical expertise leads you to deviate from those targets, you should perform testing to determine what is happening in the patient's ear canal

Can I verify using Oticon targets?

What fitting system are you using?

Which rationale do you want to use to verify?

Do you want to use AutoFit or Conventional REM?



*To perform conventional REM with VAC+/DSE, the REM module needs to be launched from the REM tool in Oticon Genie 2. Note: Generic rationales refer to non-proprietary rationales such as NAL-NL2 and DSL v5


oticon

Can I verify using Oticon targets?

What fitting system are you using?

Which rationale do you want to use to verify?

Do you want to use AutoFit or Conventional REM?



*To perform conventional REM with VAC+/DSE, the REM module needs to be launched from the REM tool in Oticon Genie 2. Note: Generic rationales refer to non-proprietary rationales such as NAL-NL2 and DSL v5

oticon

Verifying using the Speech Intelligibility Index (SII)



Assessing how much of the **speech signal is audible** to the client



Using the SII for Hearing Aid Fittings



AJA

Tutorial

The Speech Intelligibility Index: Tutorial and Applications for Children Who Are Deaf and Hard of Hearing

Kathryn B. Wiseman,^a Caitlin Sapp,^b Derek Stiles,^c Elizabeth A. Walker,^d and Ryan W. McCreery^a

^aCenter for Hearing Research, Boys Town National Research Hospital, Omaha, NE ^bDepartment of Audiology, University of North Carolina, Chapel Hill ^cDepartment of Otolaryngology and Communication Enhancement, Boston Children's Hospital, MA ^dDepartment of Communication Sciences and Disorders, University of Iowa, Iowa City

ARTICLE INFO

Article History:
Received November 25, 2024
Revision received March 27, 2025
Accepted April 14, 2025

Editor-in-Chief: Erin M. Picou
Editor: Douglas P. Sladen

https://doi.org/10.1044/2025_AJA-24-00267

ABSTRACT

Purpose: Speech audibility refers to how much access a listener has to the acoustic cues that comprise spoken communication. The Speech Intelligibility Index (SII) is a measure of speech audibility that has had increasing use in pediatric audiology to quantify auditory access with and without hearing aids. This tutorial focuses on clinical applications of the SII in children who are deaf and hard of hearing (DHH).

Method: We discuss the calculation of the SII and the current clinical tools that calculate aided and unaided SII. We describe practical uses of the SII throughout the steps of audiologic care, including at diagnosis of hearing loss, hearing technology candidacy, hearing aid verification, and long-term management of children who are DHH. Common questions about clinical applications of the SII are addressed throughout the tutorial. Finally, we present three case studies that illustrate the application of SII to clinical decision making and counseling with children and their families.

Conclusions: The SII has multiple applications in pediatric audiology. Speech audibility can support device candidacy and management decisions for children who are DHH.

audiologyonline
a Continued family site

20Q: Using the Aided Speech Intelligibility Index in Hearing Aid Fittings



Level: Intermediate | 1 Hour | Course: #31978

Course Description

This course explains the Speech Intelligibility Index (SII) and how hearing care professionals can integrate it into their hearing aid fittings. It is written in an engaging Q & A format.

Presenters

Susan Scollie, PhD

Scollie, S.D. (2018). 20Q: Using the aided speech intelligibility index in hearing aid fittings. AudiologyOnline, Article23707. Retrieved from www.audiologyonline.com

Wiseman, K. B., Sapp, C., Stiles, D., Walker, E. A., & McCreery, R. W. (2025). The Speech Intelligibility Index: Tutorial and Applications for Children Who Are Deaf and Hard of Hearing. *American journal of audiology*, 1–17. Advance online publication. https://doi.org/10.1044/2025_AJA-24-00267

70 Voice Aligned Compression (VAC+)

oticon



Evidence from a clinical investigation

Highlights from a clinical investigation

Ghamkhar, 2025.
Oticon Whitepaper

VAC+ performed comparably to NAL-NL2

25% improvement in subjective overall hearing with VAC+

Higher SII at 50dB for VAC+ than for NAL-NL2

oticon

Methods


Participants




8 first-time
hearing aid
users



Aged between
52 and **82** years
(mean = 68)




Bilateral
sensorineural
hearing loss



No known
cognitive
impairment

Methods

Trials of VAC+ and NAL-NL2



Methods

Questionnaire covering
a range of auditory aspects

Items rated on
5-point Likert scale


Open-ended
questions

Choice between
VAC+ and NAL-NL2

Results

Speech intelligibility index (SII)

NAL-NL2




Ghamkhar, 2025.
Oticon Whitepaper

oticon

Results

Speech intelligibility index (SII)



At 50 dB, **SII was higher when fit to VAC+** than NAL-NL2

Ghamkhar, 2025.
Oticon Whitepaper

oticon

Results

Questionnaires

Adaptation
time

Speech
clarity

Sound
quality

Listening
comfort

Help in **loud**
environments


No difference or preference between rationales
for the majority of **questionnaire items**

Ghamkhar, 2025.
Oticon Whitepaper

oticon

Results

Perceived improvement in overall hearing



Rationale

25%

**improvement in subjective
overall hearing** when fit to VAC+

Ghamkhar, 2025.
Oticon Whitepaper

oticon

Conclusion


VAC+ provided **better access to soft speech**

VAC+ gave a **25% improvement** in **overall hearing experience**

VAC+ is a **clinically valid rationale** and a **great alternative to NAL-NL2**

Ghamkhar, 2025.
Oticon Whitepaper

oticon



Learn more...



The Oticon VAC+ fitting rationale. [White paper]

References

- American National Standards Institute. (1997, R2020). Methods for calculation of the Speech Intelligibility Index, ANSI S3.5-1997. New York: Acoustical Society of America
- Buus, S., & Florentine, M. (2002). Growth of loudness in listeners with cochlear hearing losses: recruitment reconsidered. *Journal of the Association for Research in Otolaryngology : JARO*, 3(2), 120–139. <https://doi.org/10.1007/s101620010084>
- Byrne, D., Dillon, H., Ching, T., Katsch, R., & Keidser, G. (2001). NAL-NL1 procedure for fitting nonlinear hearing aids: characteristics and comparisons with other procedures. *Journal of the American Academy of Audiology*, 12(1), 37–51.
- Cox, R. M., & Alexander, G. C. (2002). The International Outcome Inventory for Hearing Aids (IOI-HA): psychometric properties of the English version. *International journal of audiology*, 41(1), 30–35. <https://doi.org/10.3109/14992020209101309>
- Dao, A., Folkeard, P., Baker, S., Pumford, J., & Scollie, S. (2021). Fit-to-Targets and Aided Speech Intelligibility Index Values for Hearing Aids Fitted to the DSL v5-Adult Prescription. *Journal of the American Academy of Audiology*, 32(2), 90–98. <https://doi.org/10.1055/s-0040-1718707>
- Dillon, H. (2012). *Hearing aids* (2nd ed.). Thieme.
- Ghamkhar, M. (2025). *The Oticon VAC+ fitting rationale*. [White paper]. Oticon.
- Keidser, G., Dillon, H., Flax, M., Ching, T., & Brewer, S. (2011). The NAL-NL2 Prescription Procedure. *Audiology research*, 1(1), e24. <https://doi.org/10.4081/audiores.2011.e24>
- Kitterick, P., & Edwards, B. (2025). NL3: The Next Generation Hearing Aid Fitting System. Presented at the American Academy of Audiology 2025 and HearTech Expo
- Le Goff, N. (2015). Amplifying soft sounds – a personal matter. Oticon whitepaper.
- Le Goff, N., & Schum, D. (2015). Evidence on Soft Speech Booster. Oticon whitepaper.
- Marozeau, J., & Florentine, M. (2007). Loudness growth in individual listeners with hearing losses: a review. *The Journal of the Acoustical Society of America*, 122(3), EL81. <https://doi.org/10.1121/1.2761924>

References

- National Acoustic Laboratories. (2025). NAL-NL3: The Next Generation Fitting System. https://www.nal.gov.au/nal_products/nal-nl3-the-next-generation-fitting-system/#:~:text=NAL-NL3%20is%20a%20prescriptive%20fitting%20system%20used%20by,amplification%20for%20individuals%20based%20on%20their%20hearing%20loss.
- Powers, T. A., Branda, E., & Beilin, J. (2014, June 16). Clinical comparison of a manufacturer's proprietary fitting algorithm to the NAL-NL2 prescriptive method. *AudiologyOnline*. <https://www.audiologyonline.com/articles/clinical-comparison-manufacturer-s-proprietary-12708>
- Schmidt, E. & Caporali, S. & Jensen, N.S. (2017). Match to Target: How to Get Accurate First Fits. *Hearing Review*. 24. 16-27.
- Scollie, S., Seewald, R., Cornelisse, L., Moodie, S., Bagatto, M., Lurnagaray, D., Beaulac, S., & Pumford, J. (2005). The Desired Sensation Level multistage input/output algorithm. *Trends in amplification*, 9(4), 159–197. <https://doi.org/10.1177/108471380500900403>
- Scollie, S. (2007, January 15). *DSL version v5.0: Description and Early Results in Children*. *Audiology Online*. <https://www.audiologyonline.com/articles/dsl-version-v5-0-description-959>
- Scollie, S. (2018, September 10). *20Q: Using the Aided Speech Intelligibility Index in Hearing Aid Fittings*. *Audiology Online*. <https://www.audiologyonline.com/articles/20q-aided-speech-intelligibility-index-23707>
- Spencer, K (2024). *History of the Cell Phone (1973 to 2024)*. *PracticallyNetworked.com*. <https://www.practicallynetworked.com/history-of-the-cell-phone/>
- Wiseman, K. B., Sapp, C., Stiles, D., Walker, E. A., & McCreery, R. W. (2025). The Speech Intelligibility Index: Tutorial and Applications for Children Who Are Deaf and Hard of Hearing. *American journal of audiology*, 1–17. Advance online publication. https://doi.org/10.1044/2025_AJA-24-00267

