

What is the right solution for my patient?

How to explain the differences between the automatic Program 1 in ReSound Sensia microRIE technology levels 9 and 7.

	microRIE Tech Level 9	microRIE Tech Level 7
Hearing in background noise		
Directionality 	Intelligent All-Around Provides automatic adjustments in all listening environments, including DNN noise reduction automatically activating to improve speech intelligibility and listening comfort. “The hearing aids automatically adjust to optimize to every listening environment so that you can hear what you want to hear, without needing to make manual adjustments as the environment gets more complex and noisy.”*	360 All-Around Provides natural binaural hearing with 4-mic beamforming to support speech clarity and listening effort. “The hearing aids automatically adjust to the environment to provide support in complex and noisy listening situations. A manual program change may be necessary in very complex listening environments, like a noisy restaurant.”*
Comfort in noise		
Noise Reduction 	A strong level of noise reduction for added comfort. “Uses two types of noise reduction, including smart AI technology, to automatically filter out background noise while keeping speech clear. The louder the environment, the harder it works, so listening feels easier and more comfortable.”* Up to -15dB	A considerable level of noise reduction for added comfort. “Uses a noise reduction strategy to separate speech from noise and requires a manual program change to change the amount of noise reduction.”* Up to -10 dB
Automatic adjustments		
Environmental Optimizer 	Automatic personalized volume control and noise reduction. “Adjusts the volume and noise reduction to adapt to the changes in your environment so you can be comfortable in your surroundings.”*	Automatic personalized volume control. “Adjusts the volume to adapt to the changes in your environment so you can be comfortable in your surroundings. The amount of noise reduction does not adjust automatically.”*
Personalized adjustments		
Sound Enhancer 	Adjustment tool within the ReSound Smart 3D™ app that allows patients to customize frequency bands, noise reduction and wind noise reduction for the most personalized fit. “Allows you to customize the sound to your exact preferences in any environment.”*	Adjustment tool within the ReSound Smart 3D app that allows patients to customize frequency bands only. “Allows you to customize sounds, like an equalizer on a stereo, for a more personalized sound experience.”*

*A suggestion on how to explain the product benefits to your patients.

Manufacturer according to FDA:

GN ReSound North America

5005 Dean Lakes Blvd.
Shakopee, MN 55379
USA
1-800-248-4327
pro.resound.com

Manufacturer according to Health Canada:

ReSound Canada

2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
pro.resound.com

