

Unleash the Power of Conversation


Signia Integrated Xperience (IX)
with RealTime Conversation Enhancement



signia

Group Conversation in Noise is Difficult...

- Group conversations in noise are the biggest unmet need for people with hearing Loss¹
- Can lead to frustration, isolation and ultimately withdrawal²
- They miss out on the moment, we miss out on their contribution to the conversation



¹ MarkeTrak 2022: Insights into the Evolving Hearing Market: Hearing Aid Benefit and Satisfaction Results from the MarkeTrak 2022 Survey: Importance of Features and Hearing Care Professionals - PMC (nih.gov) Picou, (2022).

² Impact of Hearing Loss on Daily Life and the Workplace - Hearing Loss - NCBI Bookshelf (nih.gov)

Signia's Integrated Xperience provides a revolutionary three-stage process that puts conversation first.

ANALYZE



Analyzes 192,000 data points every second to authentically reflect real-life conversations

AUGMENT



Leverages the comfort and clarity of split processing to provide an immersive conversational experience without distractions

ADAPT



Focus streams adapt 1,000 times per second to ensure wearers can fully engage in group conversations



Benefits of RealTime Conversation Enhancement:

A comfortable,
immersive
conversational
experience without
distractions

Access to multiple
speakers via the
world's first multi-
stream architecture¹
with speech that
pops above the
surrounding sound


Keep up and
contribute to the
conversation



95%
of wearers have
better speech
reception
with RealTime
Conversation
Enhancement
than without.²

¹ Number of focus streams varies per technology level.

² Visit signia-pro.com/en-us/claims-data for more details



Future-proof connectivity & performance



With Signia, you give your patients the chance to stay connected with the future of hearing care.

- ASHA (Audio Streaming for Hearing Aids)
- Made for iPhone
- LE Audio ready via future firmware update
- Bluetooth connectivity via streamer

Signia also provides firmware updates to enhance performance of your devices.

Own Voice Processing 2.0



Own Voice Processing (OVP) promotes hearing wellness

Contribute comfortably — OVP ensures own voice comfort, activating whenever you speak to keep your own voice sounding natural and clear.


Boost Communication — 78% of participants reported increased communication when using OVP.³

Signia Assistant



The 24/7 hearing companion

Signia Assistant is the hearing aid wearer's **AI-powered companion**, understanding your environment and adapting the hearing aids for personalized support.⁴



³ Høydal EH. A new own voice processing system for optimizing communication. Hearing Review.

⁴ Accessible via My WellBeing in the Signia App.

Pure Charge&Go IX / T IX

Deliver discreet, advanced solutions

Pure Charge&Go IX and Pure Charge&Go T IX offer advanced hearing technology in an attractive and comfortable receiver-in-canal (RIC) design.

Power through the day

A rechargeable battery ensures your clients enjoy all-day hearing support.
Battery run time up to 39 hours⁵

Stay connected

- ASHA (Audio Streaming for Hearing Aids)
- Made for iPhone®
- LE Audio ready
- Bluetooth® connectivity via streamer
- Pure Charge&Go T IX comes with a T-coil

Proven reliability ensuring best performance

- Proven water and oil protection
- IP68 rating ensures the hearing aids will still function if immersed under 1 m of water for 1h⁶


Available Chargers:



Pure Portable
Charger



Pure
Charger



Pure Dry&Clean
Charger

Available Colors:



Black



Graphite



Dark
Champagne



Silver



Pearl
White



Fine
Gold



Deep
Brown



Sandy
Brown



Rose
Gold



Beige

⁵ Pure C&G T IX battery run time up to 39 hours without streaming or 34 hours with 5 hours streaming. Pure C&G IX battery run time up to 24 hours including 5 hours streaming.

⁶ Excluding receiver

Silk Charge&Go IX

Seamless hearing with discretion

Silk Charge&Go IX – an innovative completely-in-canal (CIC) hearing aid that pairs exceptional sound quality with unmatched discretion.

Power through the day

The convenience of a rechargeable battery for all-day listening. Plus, the accompanying hassle-free case charges wirelessly.

Battery run time up to 28 hours

Conveniently charge hearing aids up to 4 times on the go.

Binaural OneMic Directionality

The only CIC device that offers Binaural OneMic Directionality to focus on a speaker in noisy listening environments.

Unmatched discretion

The soft, textured surface avoids reflections for added discretion.

Immediate fit, immediate impact

The quick, satisfying experience of a ready-to-wear CIC. The Silk Charge&Go IX's small, snug design is discreet and it can be worn immediately.

Control with the app

The user-friendly Signia app and miniPocket remote control allows patients to fine-tune their hearing aids to their unique needs for personalized control anytime, anywhere.

Available Faceplate Colors:



Black



Mocha



Learn more today.

Made for
 iPhone

"Made for iPhone" mean that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone may affect wireless performance. iPad is a trademark of Apple Inc., registered in the U.S. and other countries.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc.

signia.net | signia-pro.com

Signia is a registered trademark of WSAUD A/S. SIV-17175

Be
Brilliant™