

Multi-Mic+

The world's first remote mic with Auracast™

The new Multi-Mic+*, enabled with Auracast broadcast audio, helps cut through the noise and keeps your conversations going when situations get extra loud. The audio picked up by the Multi-Mic+ is streamed directly to your ReSound Nexia hearing aids giving you extra support in situations where you struggle to hear.



Presenter mode

One-to-one or one-to-many, clip the Multi-Mic+ vertically onto the clothing of your conversation partner or a presenter when you're at a noisy venue, or the person who is giving a speech or presentation from across a large room.



Table mode

To hear group conversations, place the Multi-Mic+ horizontally on a table when many people are talking and you want to hear all their voices clearly, like a work meeting or a family dinner.



Audio mode

Stream audio from your personal electronic devices, like laptops, TVs and more using the headphone mode; Connect in public places where hearing loops or FM systems are available.

The world's first remote mic with Auracast

The Multi-Mic+ offers the new standard of Bluetooth connectivity. With Auracast, stream and share your audio with a larger group. Educators, presenters, event hosts and tour guides can wear the Multi-Mic+ and broadcast their sound experience to other compatible listeners.

Auracast allows the Multi-Mic+ to broadcast audio simultaneously to an unlimited number of nearby compatible audio devices with Bluetooth® Low Energy (LE) Audio. In comparison, traditional microphone technology can only stream audio to hearing instruments of the same brand.



Eliminate noise with the new Multi-Mic+



Multi-Mic+ Features

- Small, sleek and fingerprint resistance design with optimized buttons for easy handling
- Long battery life, offering 15 hours in mic mode and up to 20 hours when streaming audio directly from an audio source like a computer
- 82 feet (25 meters) of streaming range
- Easy OneTouch pairing by simply pressing the pairing button on the back of the Multi-Mic+ and placing the hearing aids within 2 inches (5 centimeters) of the Multi-Mic+
- Share exciting audio experiences in both private and public spaces with Auracast



ReSound Smart 3D app

Conveniently adjust the Multi-Mic+ volume settings and Auracast broadcast audio transmitter settings in the ReSound Smart 3D™ app.



For more information, gs.resound.com

Manufacturer according to FDA:

ReSound Government Services

8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
gs.resound.com



View instructional videos on
youtube.com/resoundus.