

CROS Transmitter

A solution for patients with single-sided deafness



A complete hearing solution

Bring the Resound Nexia™ experience to even more patients with the wireless CROS/BiCROS solution - now compatible with RIEs, BTEs and rechargeable customs.

- Available in the industry's smallest wireless Micro RIE CROS and our new wireless BTE CROS
- Provides all-day* battery life on one charge while streaming with CROS functionality
- Intuitive CROS balance control within ReSound Smart Fit™
- Convenient patient access to CROS balance control within ReSound Smart 3D™ app
- Connects with Bluetooth® Low Energy (LE) Audio and Auracast™
- Compatible with all wireless accessories including TV-streamer+ and Multi-Mic+

ReSound Nexia RIE CROS



CX160S-DRWC

compatible with RIE models


ReSound Nexia BTE CROS



CX177-DWC

compatible with rechargeable customs and BTE models

Available colors




*All-day battery life on one charge is based on typical use time of 16 hours (both transmitter and hearing aid). Battery life can vary by use.

Fully integrated in the ReSound Smart Fit 2.0 fitting software

The solution seamlessly integrates into the fitting flow, allowing for easy adjustments to CROS/BiCROS fittings.

Additionally, the fitting software provides the flexibility to assign preferred controls to individual patients, ensuring a personalized and intuitive experience.



ReSound Smart 3D app

ReSound Nexia CROS is a complete solution offering direct audio streaming and more through the ReSound Smart 3D™ app. Patients can conveniently adjust CROS balance control within the app at any time.



Compatibility & Features

	Micro RIE CROS	BTE CROS
Compatibility	All ReSound Nexia RIEs	ReSound Nexia BTEs and rechargeable customs
Battery life	16 hours*	16 hours*
Charger options	Premium, desktop (CFF 8)	Desktop (CFF 70)
Fitting options	SureFit 3 LP, MP, HP receivers with SureFit 3 domes or earmold	SureFit 3 thin tubes or earhook with SureFit 3 domes or earmold
Telecoil	N/A	Yes
TapControl	Yes	N/A

*All-day battery life on one charge is based on typical use time of 16 hours (both transmitter and hearing aid). Battery life can vary by use.

Manufacturer according to FDA:

ReSound Government Services
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
gs.resound.com

Learn more at gs.resound.com