

Available colors

ReSound Nexia customs



ReSound custom shell only colors




ReSound Vivia (RIE)
ReSound Enzo IA (SPBTE)



*Known as Champagne in ReSound Nexia

ReSound Nexia (RIE, BTE)



ReSound ENZO Q (BTE)



Apps



Chargers



Wireless accessories



*Multi-Mic+ and TV-Streamer+ are compatible with ReSound Vivia, ReSound Nexia and ReSound Enzo IA.



Quick reference guide

Government Services
November 1, 2025 - April 30, 2026




ReSound Vivia™
ReSound Nexia™
ReSound Enzo™ IA
ReSound ENZO Q™












Quick reference support
ReSound Government Services: 800-392-9932
Government Services Consumer Support Line: 888-879-7036



Manufacturer according to FDA:
ReSound Government Services
5005 Dean Lakes Blvd.
Shakopee, MN 55379
USA
1-800-392-9932
gs.resound.com



© 2025 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc.

ReSound Vivia™		
Models	Fitting range	Max Output/Gain
 microRIE 60S Lithium-ion		M&RIE: 115/61 Low power: 113/56 Medium power: 115/60 High power: 120/60 Ultra power: 129/71
 Standard RIE 61 312		
 Standard RIE 62 13		
 microRIE 60S CROS Transmitter Lithium-ion	Compatible with: Lower power Medium power High power ReSound Vivia microRIE CROS transmitter is only compatible with ReSound Vivia RIE styles.	

ReSound Nexia™		
Models	Fitting range	Max Output/Gain
 microRIE 60S Lithium-ion		M&RIE: 115/53 Low power: 114/54 Medium power: 116/59 High power: 120/63 Ultra power: 130/73
 Standard RIE 61 312		
 Standard RIE 62 13		
 microRIE 60S CROS Transmitter Lithium-ion	Compatible with: Lower power Medium power High power ReSound Nexia microRIE CROS transmitter is only compatible with ReSound Nexia RIE styles.	
 Standard BTE 77 Lithium-ion		Open: 127/53 Closed: 128/60
 Power BTE 88* Plastic or Metal Hook Lithium-ion		Power: 133/67 High power: 133/73 *BTE 88 can be configured with plastic or metal hook in ReSound Smart Fit™

ReSound Nexia™		
Models	Fitting range	Max Output/Gain
 Standard BTE 77 CROS Transmitter Lithium-ion	Compatible with: Plastic hook Thin-Tube ReSound Nexia BTE CROS transmitter is only compatible with ReSound Nexia BTE and rechargeable custom styles.	
 Half Shell/ITC** Lithium-ion		Low power: 114/40** Medium power: 118/50 High power: 121/60 Ultra power: 130/70 **Devices will be built as small as possible based on the patient's ear anatomy. Low power available for ITC orders only.
 Full Shell Lithium-ion		
 CIC*** 10A		Low power: 114/40 Medium power: 116/50 High power: 121/60 ***CIC is non-wireless.
 ITC 312		Low power: 114/40 Medium power: 117/50 High power: 120/60 Ultra power: 130/68

ReSound Enzo™ IA • Technology Levels 9, 7, 5, 4		
Models	Fitting range	Max Output/Gain
 Super Power BTE 98 Metal Hook Lithium-ion		140/84
 Super Power BTE 98 Plastic Hook Lithium-ion		140/79

ReSound ENZO Q™ • Technology Levels 9, 7, 5		
Models	Fitting range	Max Output/Gain
 High Power BTE 88 Metal Hook 13		134/73
 Super Power BTE 98 675		141/83

Fitting Range Key
 LP and MP: Mild - Moderate
 HP: Moderate - Severe
 UP: Severe - Profound

 = Telecoil default

Max Output/Gain: 2cc Coupler