



Cochlear® Osia® 3 Sound Processor

Hear the difference

Help your patients feel unstoppable in every moment

Built on the latest Osia® technology, the Osia 3 Sound Processor elevates hearing performance with a larger dynamic range and intelligent features designed for real-world listening—empowering patients to hear more clearly, comfortably, and confidently.

Proven outcomes for your patients

The new ForwardFocus feature directs the microphones toward sounds originating in front of the listener, delivering up to a **7% improvement in speech understanding in noisy environments**.¹

The Osia 3 Sound Processor has a fitting range of **up to 65 dB HL**.^{2,3*}

Participants report **significantly less hearing difficulty** with the Osia 3 Sound Processor compared with the Osia 2 Sound Processor utilizing SSQ-12 ($p < 0.025$).¹

Based on clinical data, recipients with the Osia 3 Sound Processor have on average **4 dB more available gain** compared to the Osia 2 Sound Processor.⁴



+7%

Improvement in noise

≤65
dB HL


Fitting range

+4
dB


Available gain

Participants reported higher satisfaction when using the Osia 3 Sound Processor compared to the Osia 2 Sound Processor¹


Overall sound quality





Hearing in background noise



Listening to music



Hearing in wind



Provide your patients with the latest technology

Fully rechargeable⁵

- The rechargeable Osia 3 Sound Processor can reduce disposable battery use by approximately 200 batteries per year, compared to the Osia 2 Sound Processor
- The Osia 3 Sound Processor offers significant battery-related cost savings, estimated at \$150 USD per year and \$750 USD over 5 years

Broadest range of streaming options

- Hands-free calling and streaming on Android™ and iOS
- Bluetooth® LE Audio and Auracast[^]

Waterproofed

- The Osia 3 Sound Processor is resistant to dust and water and fulfills an IP68 waterproof down to 1 meter for 1 hour ingress protection classification[†]

We're here to help

We can help your patients navigate the insurance process. If your patients are concerned about cost, we can help them understand their options, which may include financing the out-of-pocket cost into no-interest monthly payments. Please encourage your patients to contact us at **1-800-483-3123** or email upgrade@cochlear.com for help with how to pay for their sound processor replacement.

1. Hearing Outcomes With the Osia 3 Sound Processor Compared to the Osia 2 Sound Processor in Adult Osia Implant Recipients (SONUS). Clinicaltrials.gov NCT04754477. Available at: <https://clinicaltrials.gov/study/NCT07156461?term=Osia%203%20Sound%20Processor&rank=2>
2. Cochlear Limited. D2304775-V1 Osia 3 SP and OSI200 PROF Datasheet. 2026.
3. Cochlear Limited. D2304776-V1 Osia 3(I) SP and OSI300 PROF Datasheet. 2026.
4. Cochlear Limited. D1821490. Technical Report on Available Gain - CBAS Portfolio. 2026.
5. D2423880, Osia 3 and 3(I) Sound Processors Battery Cost of Ownership mode
6. Cochlear Limited. D2304775-V1_EN Osia 3 SP and OSI200 PROF Datasheet
7. Cochlear Limited. D2304776-V1_EN Osia 3(I) SP and OSI300 PROF Datasheet

* Refer to the relevant Osia System datasheet for additional fitting range information.

[^] For compatibility information regarding Cochlear products, apps and services with certain devices and operating systems, visit www.cochlear.com/compatibility.

~ LE Audio is not compatible with certain bilateral and bimodal configurations. Please refer to the compatibility details here: <https://cochlear.com/leaudio-compatibility>

[†] The Osia 3 Sound Processor meet the IP68 rating of the International Standard IEC60529 of freshwater waterproof. These processor configurations were tested by continuous submersion in freshwater for 60 minutes at a depth of 1 meter and functioned as intended. Cochlear offers the Aqua+ accessory for additional protection during extended water use, in salty or rushing water environments. For additional information, please refer to the appropriate user guide.^{6,7}

This material is intended for health professionals. If you are a consumer, please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always read the instructions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information.

©Cochlear Limited 2026. All rights reserved. Cochlear, Hear now. And always, Nucleus, Kanso, Baha, Osia, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of the Cochlear group of companies (unless otherwise noted).

