



Cochlear®

Hear now. And always

**Cochlear® Osia® 3
Sound Processor**

**Feel unstoppable
in every moment**

**Now fully
rechargeable**

Feel unstoppable in every moment

We believe sound is more than a sense. It's the force that drives us forward: empowered, bold and full of potential.

Everyone's hearing story is unique. From conferences to cafés, favorite songs to family celebrations, every moment sounds different. The Cochlear® Osia® 3 Sound Processor introduces advanced technology to help you hear more clearly, reduce background noise, and stay connected to the voices that matter most, so you can feel confident and unstoppable in every moment.



Hearing never looked better

With the new Osia 3 Sound Processor you can express your style your way, with six new bright color covers you can switch between whenever you like.



“[The Osia 3 Sound Processor] has improved my social experiences with people. It gives me more confidence.”

- Recipient in Osia 3 Sound Processor clinical study¹



Clear sound. Bright future

The Osia 3 Sound Processor's excellent sound quality^{1,2} is driven by improvements to our unique sound technology. Its high power and clarity deliver even better sound quality than ever before.^{1,2}

You can pick up the details that matter, cut through distractions, and focus on what you do want to hear. The Osia System handles the sound, so you can enjoy the ride and live every moment with confidence.

Tuned for speech understanding

The Osia System is the most powerful device in its category³⁻⁵, driven by its unique Piezo Power™ technology and precision sound processing. High power with a large dynamic range of sound helps speech sound clearer, even in noise.^{6,7}




ForwardFocus

The new ForwardFocus feature directs microphones to focus on what's in front of you, for higher word recognition in noisy situations.¹



Excellent hearing experience

In a study, 9 of 10 recipients enjoyed or highly enjoyed their overall hearing experience in real-world situations with the Osia 3 Sound Processor.¹



Large dynamic range

In a hearing device, a large "dynamic range" (the span between the softest and loudest sounds) means soft sounds are amplified while louder sounds remain crisp and less distorted.^{6,7}



Better hearing outcomes

In a study, recipients with the Osia 3 Sound Processor showed significantly improved hearing performance compared to the previous generation.¹



**"The noisy background faded,
and I could focus on
what I needed to."**

- Recipient in Osia 3 Sound Processor clinical study¹

Hear the difference.

In a clinical study, participants reported a higher satisfaction when using the Osia 3 Sound Processor compared to the Osia 2 Sound Processor.

Results:



Scan the QR code to see a more detailed comparison on how the Osia 3 Sound Processor compares with the previous generation.

Splash-proof. Ready for your life

With the Osia 3 Sound Processor, you can say “Yes!” to that spa day and lakeside vacation. With **best-in-class dust and water protection** at level IP68[‡], you can do more of what you love.



Henry, Cochlear
Osia System recipient



Tested and certified

As the only IP68-certified device in its category³⁻⁵, the Osia 3 Sound Processor is protected for everyday sweat, unexpected rains, a dip in the bathtub or a short swim in freshwater.[‡]



Dive deeper

Add Aqua+ to swim, snorkel or surf with confidence, knowing your sound processor is protected for saltwater and longer underwater exposure.[‡] And for extra peace of mind, choose from a range of optional retention accessories.



One charge. All-day confidence

Experience sound that moves with your day. As the first fully rechargeable sound processor in its class, the Osia 3 Sound Processor keeps you connected from early mornings to late-night moments.



Up to 30 hours of hearing

Hear your best all day and night with up to 30 hours on a single charge.†



On-the-go charging case

Similar to wireless earbuds, it comes with a charging case that holds up to two full charges. So you'll also be covered for that well-deserved weekend getaway.



Quick charging

When your busy life gets ahead of you, just 15 minutes of charging can get you up to two hours of hearing. Fully charge both the sound processor and the case in three hours.



No more batteries

No need to buy, store and keep track of tiny disposable batteries that lead to more hassle and household waste.

Powerful performance that fits
modern expectations.

Effortless streaming. Endless connections

With your Osia 3 Sound Processor, you can step into a more connected world. With the broadest range of streaming options available in a bone conduction device⁺, it's designed to keep you connected to what matters wherever life takes you.

iOS and Android™ streaming

Stay effortlessly connected to the songs, stories and moments you love.

Bluetooth® LE Audio enabled

Stream high-quality sound from compatible phones, laptops and other devices with the most advanced Bluetooth platform.^{8,9}

Auracast™ broadcast audio enabled

Share music with a friend and connect to public audio streams, as more devices and venues adopt Auracast broadcast audio.

Hands-free calling

Use your sound processor as a headset and take calls in quiet environments without holding up your phone to speak.^Δ

Learn more at www.cochlear.com/Connectivity



NEW



NEW



Osia Smart App⁺



Stay in control: Easily and discreetly adjust your sound processor settings, change programs and activate ForwardFocus with the Osia Smart App.



Streaming Hub: Easily manage your wireless devices, including the new Multi-Mic+ and TV-Streamer+ with one-touch pairing and automatic streaming.



Remote Care: Connect with your clinician from home via a Remote Assist video appointment, so they can fully program and adjust your sound processor: no travel needed.



Hearing Tracker: Automatically track your daily sound processor usage, so you and your audiologist can make decisions about your care that are backed by data.

Hear now. And always

Cochlear is dedicated to helping people with moderate to profound hearing loss experience a world full of hearing. As the global leader in implantable hearing solutions, we have helped more than 750,000 people of all ages to hear and connect with life's opportunities.

We aim to give people the best lifelong hearing experience and access to next generation technologies. We collaborate with leading clinical, research and support networks to advance hearing science and improve care.

That's why more people choose Cochlear than any other hearing implant company.

References

1. CLTD5861 Clinical Investigation Report. Cochlear Limited, Sweden. 2026; D2367205.
2. CLTD5839 Clinical Investigation Report. Cochlear Limited, Sweden. 2025; D2367209.
3. Osia 3(I) Sound Processor and OSI300 Implant Datasheet. Cochlear Limited, Sweden. 2026; D2304776.
4. BCI602 BONEBRIDGE System Factsheet. Med-El. 29214 r4.0.
5. Sentio 1 Mini Product Information. Oticon Medical. 275144en – Version A / 2024.04.
6. Van Eeckhoutte M, Folkeard P, Glista D, Scollie S. Speech recognition, loudness, and preference with extended bandwidth hearing aids for adult hearing aid users. *Int J Audiol.* 2020 20;1-12.
7. Gawliczek T, Wimmer W, Caversaccio M, Kompis M. Influence of maximum power output on speech understanding with bone anchored hearing systems. *Acta Otolaryngol.* 2020;140(3):225-229.
8. Hunn, N. Introducing Bluetooth® LE Audio. Version 2. November 2024. [Cited July 2025]. Available from: <https://www.bluetooth.com/bluetooth-resources/le-audio-book-v2>.
9. Afaneh, M. A technical overview of LC3 [Cited July 2025] <https://www.bluetooth.com/blog/a-technical-overview-of-lc3>.
10. Cochlear Limited. D2304775-V1_EN Osia 3 SP and OSI200 PROF Datasheet.
11. Cochlear Limited. D2304776-V1_EN Osia 3(I) SP and OSI300 PROF Datasheet.

† Battery life varies for every user, according to the age of the battery, the programs used each day, your implant type, the thickness of skin covering your implant, and the size and type of battery used. Streaming may decrease sound processor battery life depending on how often and for how long streaming is engaged.

+ For information regarding the compatibility of Cochlear's sound processors and Smart Apps with Apple devices, Android devices, or our Remote Care services, visit www.cochlear.com/compatibility

^ LE Audio is not compatible with certain bilateral and bimodal configurations. Please refer to the compatibility details here: <https://cochlear.com/leaudio-compatibility>

Δ Hands-free calling is intended for use in quiet environments. See User Manual for more details.

‡ The Osia 3 Sound Processor meet the IP68 rating of the International Standard IEC60529 of freshwater waterproof. These processor configurations were tested by continuous submersion in freshwater for 60 minutes at a depth of 1 meter and functioned as intended. Cochlear offers the Aqua+ accessory for additional protection during extended water use, in salty or rushing water environments. For additional information, please refer to the appropriate user guide.^{10,11}

The Cochlear Osia Smart App is available on App Store and Google Play. For compatibility information visit www.cochlear.com/compatibility

The Cochlear Osia 3 and 3(I) Sound Processors are compatible with Apple devices. For compatibility information visit www.cochlear.com/compatibility.

Please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always read the instructions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information.

Views expressed are those of the individual. Consult your hearing health provider to determine if you are a candidate for Cochlear technology. Outcomes and results may vary.

Android is a trademark of Google LLC. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Google Play and the Google Play logo are trademarks of Google LLC.

Apple, the Apple logo, Apple Watch, FaceTime, Made for iPad logo, Made for iPhone logo, Made for iPod logo, iPhone, iPad Pro, iPad Air, iPad mini, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear is under license.

As Bluetooth LE Audio compatible devices become available, a sound processor firmware update will be required to use certain features.

Auracast™ broadcast audio capability is subject to third party adoption of the Auracast protocol.

The Bluetooth® and Auracast™ word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear is under license.

©Cochlear Limited 2026. All rights reserved. ACE, Advance OffStylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントゥア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, Nexa, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, SoundBand, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield and Xidium are either trademarks or registered trademarks of the Cochlear group of companies.

Cochlear North America

10350 Park Meadows Drive, Lone Tree, CO 80124, USA
Tel: +1 303 790 9010 Support: Tel: +1 800 483 3123

Cochlear Canada Inc.

2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada
Tel: +1 800 483 3123 Fax: +1 416 972 5083

www.cochlear.com



OSI261 ISS1 JUL26

Follow us     