

ReSound Sensia™ CROS Transmitter


A solution for patients
with single-sided deafness



A complete hearing solution

Bring the ReSound Sensia experience to even more patients with the new wireless CROS/ BiCROS solution. Available in the microRIE style and color options, the ReSound Sensia CROS Transmitter perfectly pairs with all ReSound Sensia RIE models. It delivers premium CROS and BiCROS performance, designed for exceptional comfort and all-day confidence, with the lowest internal noise¹ for more effortless listening.

- Available in the industry's smallest rechargeable CROS² design
- Compatible with ReSound Sensia RIE models
- Programmable as a CROS or BiCROS solution
- Provides all-day* battery life on one charge while streaming with CROS functionality
- Intuitive CROS balance control within ReSound Smart Fit™
- Convenient patient access to CROS balance control within ReSound Smart 3D™ app
- Connects with Bluetooth® Low Energy (LE) Audio and Auracast™
- Compatible with all wireless accessories including TV-streamer+ and Multi-Mic+



CROS Transmitter
SN060S-DRWC

Compatible styles



microRIE SNx60S-DRWC
RIE 61 SNx61-DRW
RIE 62 SNx62-DRW

Colors



¹Schumacher & Groth (2026)

²GN Proprietary data on file (2025)

*All-day battery life on one charge is based on typical use time of 16 hours (both transmitter and hearing aid). Battery life can vary by use.

Fully integrated in ReSound Smart Fit fitting software

The solution seamlessly integrates into the fitting flow, allowing for easy adjustments to CROS/BiCROS fittings.

Additionally, the fitting software provides the flexibility to assign preferred controls to individual patients, ensuring a personalized and intuitive experience.

CROS Balance

Slide to adjust the sound balance between the CROS and hearing aid

Hearing aid microphone

Activate or deactivate the sound from the hearing aid side

CROS indicator

On left in example shown

Mute

Silence transmitting signal from CROS mic

CROS indicator

On left in example shown

Mute

Silence transmitting signal from CROS mic

ReSound Smart 3D app

Patients can enjoy the latest Bluetooth® LE Audio, direct audio streaming, Auracast™ access, and intuitive everyday control through the ReSound Smart 3D app.

Download on the App Store

GET IT ON Google Play

Compatibility & Features

ReSound Sensia microRIE CROS	
Compatibility	All ReSound Sensia RIE models
Battery life	16 hours*
Charger options	Premium, standard, desktop** (CFF 8)
Fitting options	SureFit 3 LP, MP, HP receivers with domes or earmold
Telecoil	N/A
TapControl	Yes

*All-day battery life on one charge is based on typical use time of 16 hours (both transmitter and hearing aid). Battery life can vary by use.
**Desktop charger is only available in the US for RIE styles.



Manufacturer according to FDA:

GN ReSound North America
5005 Dean Lakes Blvd.
Shakopee, MN 55379
USA
1-800-248-4327
pro.resound.com

Manufacturer according to Health Canada:

ReSound Canada
2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
pro.resound.com

Learn more at pro.resound.com